

according to Regulation (EC) No. 1907/2006 as amended by (EC) No. 1272/2008

### Section 1. Identification of the Substance/Mixture and of the Company/Undertaking

- 1.1 Product Code:** 20935  
**Product Name:** 5-bromo-5-nitro-1,3-Dioxane  
**Synonyms:** Bronidox;
- 1.2 Relevant identified uses of the substance or mixture and uses advised against:**  
**Relevant identified uses:** For research use only, not for human or veterinary use.
- 1.3 Details of the Supplier of the Safety Data Sheet:**
- Company Name:** Cayman Chemical Company  
1180 E. Ellsworth Rd.  
Ann Arbor, MI 48108
- Web site address:** www.caymanchem.com
- Information:** Cayman Chemical Company +1 (734)971-3335
- 1.4 Emergency telephone number:**
- Emergency Contact:** CHEMTREC Within USA and Canada: +1 (800)424-9300  
CHEMTREC Outside USA and Canada: +1 (703)527-3887

### Section 2. Hazards Identification

- 2.1 Classification of the Substance or Mixture:**  
**Acute Toxicity: Oral, Category 4**  
**Skin Corrosion/Irritation, Category 1A-1C**  
**Serious Eye Damage/Eye Irritation, Category 1**  
**Skin Sensitization, Category 1**  
**Aquatic Toxicity (Chronic), Category 1**

**2.2 Label Elements:**



**GHS Signal Word:** **Danger**

**GHS Hazard Phrases:**

H302: Harmful if swallowed.  
H314: Causes severe skin burns and eye damage.  
H317: May cause an allergic skin reaction.  
H318: Causes serious eye damage.  
H410: Very toxic to aquatic life with long lasting effects.

**GHS Precaution Phrases:**

P260: Do not breathe {dust/fume/gas/mist/vapors/spray}.  
P264: Wash {hands} thoroughly after handling.  
P272: Contaminated work clothing should not be allowed out of the workplace.  
P273: Avoid release to the environment.  
P280: Wear {protective gloves/protective clothing/eye protection/face protection}.

**GHS Response Phrases:**

P301+312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.  
P302+352: IF ON SKIN: Wash with plenty of soap and water.  
P303+361+353: IF ON SKIN (or hair): Remove/take off immediately all contaminated clothing. Rinse skin with

water/shower.

P304+340: IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310: Immediately call a POISON CENTER or doctor/physician.

P321: Specific treatment {see ... on this label}.

P330: Rinse mouth.

P333+313: If skin irritation or rash occurs, seek medical advice/attention.

P362+364: Take off contaminated clothing and wash it before reuse.

P363: Wash contaminated clothing before reuse.

P391: Collect spillage.

### GHS Storage and Disposal Phrases:

Please refer to Section 7 for Storage and Section 13 for Disposal information.

**2.3 Adverse Human Health** Causes severe skin burns and serious eye damage.

**Effects and Symptoms:** Harmful if swallowed.

Material may be irritating to the mucous membranes and upper respiratory tract.

May be harmful by inhalation or skin absorption.

May cause an allergic skin reaction.

May cause eye, skin, or respiratory system irritation.

Very toxic to aquatic life with long lasting effects.

To the best of our knowledge, the toxicological properties have not been thoroughly investigated.

## Section 3. Composition/Information on Ingredients

CAS # / RTECS #	Hazardous Components (Chemical Name)/ REACH Registration No.	Concentration	EC No./ EC Index No.	GHS Classification
30007-47-7 JG9650000	Bromonitrodioxane	100.0 %	250-001-7 NA	Acute Tox.(O) 4: H302 Skin Corr. 1A-C: H314 Skin Sens. 1: H317 Aquatic (C) 1: H410 Eye Damage 1: H318

## Section 4. First Aid Measures

### 4.1 Description of First Aid Measures:

**In Case of Inhalation:** Remove to fresh air. If not breathing, give artificial respiration or give oxygen by trained personnel. Get immediate medical attention.

**In Case of Skin Contact:** Immediately wash skin with soap and plenty of water for at least 15 minutes. Remove contaminated clothing. Get medical attention if symptoms occur. Wash clothing before reuse.

**In Case of Eye Contact:** Hold eyelids apart and flush eyes with plenty of water for at least 15 minutes. Have eyes examined and tested by medical personnel.

**In Case of Ingestion:** Wash out mouth with water provided person is conscious. Never give anything by mouth to an unconscious person. Get medical attention. Do NOT induce vomiting unless directed to do so by medical personnel.

### Section 5. Fire Fighting Measures

- 5.1 Suitable Extinguishing Media:** Use alcohol-resistant foam, carbon dioxide, water, or dry chemical spray.  
**Media:** Use water spray to cool fire-exposed containers.  
**Unsuitable Extinguishing Media:** A solid water stream may be inefficient.
- 5.2 Flammable Properties and Hazards:** No data available.  
**Hazards:** No data available.  
**Flash Pt:** No data.  
**Explosive Limits:** LEL: No data. UEL: No data.  
**Autoignition Pt:** No data.
- 5.3 Fire Fighting Instructions:** As in any fire, wear self-contained breathing apparatus pressure-demand (NIOSH approved or equivalent), and full protective gear to prevent contact with skin and eyes.

### Section 6. Accidental Release Measures

- 6.1 Protective Precautions, Protective Equipment and Emergency Procedures:** Avoid raising and breathing dust, and provide adequate ventilation.  
As conditions warrant, wear a NIOSH approved self-contained breathing apparatus, or respirator, and appropriate personal protection (rubber boots, safety goggles, and heavy rubber gloves).
- 6.2 Environmental Precautions:** Take steps to avoid release into the environment, if safe to do so.
- 6.3 Methods and Material For Containment and Cleaning Up:** Contain spill and collect, as appropriate.  
Transfer to a chemical waste container for disposal in accordance with local regulations.

### Section 7. Handling and Storage

- 7.1 Precautions To Be Taken in Handling:** Avoid breathing dust/fume/gas/mist/vapours/spray.  
Avoid prolonged or repeated exposure.
- 7.2 Precautions To Be Taken in Storing:** Keep container tightly closed.  
Store in accordance with information listed on the product insert.

### Section 8. Exposure Controls/Personal Protection

- 8.1 Exposure Parameters:**
- 8.2 Exposure Controls:**
- 8.2.1 Engineering Controls (Ventilation etc.):** Use process enclosures, local exhaust ventilation, or other engineering controls to control airborne levels below recommended exposure limits.
- 8.2.2 Personal protection equipment:**
- Eye Protection:** Safety glasses
- Protective Gloves:** Compatible chemical-resistant gloves
- Other Protective Clothing:** Lab coat
- Respiratory Equipment (Specify Type):** NIOSH approved respirator, as conditions warrant.
- Work/Hygienic/Maintenance Practices:** Do not take internally.  
Facilities storing or utilizing this material should be equipped with an eyewash and a safety shower.  
Wash thoroughly after handling.  
No data available.

### Section 9. Physical and Chemical Properties

#### 9.1 Information on Basic Physical and Chemical Properties

<b>Physical States:</b>	<input type="checkbox"/> Gas	<input type="checkbox"/> Liquid	<input checked="" type="checkbox"/> Solid
<b>Appearance and Odor:</b>	A solid		
<b>pH:</b>	No data.		
<b>Melting Point:</b>	No data.		
<b>Boiling Point:</b>	No data.		
<b>Flash Pt:</b>	No data.		
<b>Evaporation Rate:</b>	No data.		
<b>Flammability (solid, gas):</b>	No data available.		
<b>Explosive Limits:</b>	LEL: No data.		UEL: No data.
<b>Vapor Pressure (vs. Air or mm Hg):</b>	No data.		
<b>Vapor Density (vs. Air = 1):</b>	No data.		
<b>Specific Gravity (Water = 1):</b>	No data.		
<b>Solubility in Water:</b>	No data.		
<b>Octanol/Water Partition Coefficient:</b>	No data.		
<b>Autoignition Pt:</b>	No data.		
<b>Decomposition Temperature:</b>	No data.		
<b>Viscosity:</b>	No data.		

#### 9.2 Other Information

<b>Percent Volatile:</b>	No data.		
<b>Molecular Formula &amp; Weight:</b>	C4H6BrNO4	212.0	

### Section 10. Stability and Reactivity

<b>10.1 Reactivity:</b>	No data available.		
<b>10.2 Stability:</b>	Unstable <input type="checkbox"/>	Stable <input checked="" type="checkbox"/>	
<b>10.3 Stability Note(s):</b>	Stable if stored in accordance with information listed on the product insert.		
<b>Polymerization:</b>	Will occur <input type="checkbox"/>	Will not occur <input checked="" type="checkbox"/>	
<b>10.4 Conditions To Avoid:</b>	No data available.		
<b>10.5 Incompatibility - Materials To Avoid:</b>	strong oxidizing agents		
<b>10.6 Hazardous Decomposition or Byproducts:</b>	carbon dioxide carbon monoxide hydrogen bromide gas nitrogen oxides		

### Section 11. Toxicological Information

<b>11.1 Information on Toxicological Effects:</b>	The toxicological effects of this product have not been thoroughly studied. 5-bromo-5-nitro-1,3-Dioxane - Toxicity Data: Oral LD50 (rat): 455 mg/kg; Intraperitoneal LD50 (rat): 31 mg/kg; Oral LD50 (mouse): 590 mg/kg;
<b>Chronic Toxicological Effects:</b>	5-bromo-5-nitro-1,3-Dioxane - Investigated as a primary irritant. Only select Registry of Toxic Effects of Chemical Substances (RTECS) data is presented here. See actual entry in RTECS for complete information. 5-bromo-5-nitro-1,3-Dioxane RTECS Number: JG9650000

CAS #	Hazardous Components (Chemical Name)	NTP	IARC	ACGIH	OSHA
30007-47-7	Bromonitrodioxane	n.a.	n.a.	n.a.	n.a.

### Section 12. Ecological Information

<b>12.1 Toxicity:</b>	Avoid release into the environment. Runoff from fire control or dilution water may cause pollution.
<b>12.2 Persistence and Degradability:</b>	No data available.
<b>12.3 Bioaccumulative Potential:</b>	No data available.
<b>12.4 Mobility in Soil:</b>	No data available.
<b>12.5 Results of PBT and vPvB assessment:</b>	No data available.
<b>12.6 Other adverse effects:</b>	No data available.

### Section 13. Disposal Considerations

<b>13.1 Waste Disposal Method:</b>	Dispose in accordance with local, state, and federal regulations.
------------------------------------	---

### Section 14. Transport Information

<b>14.1 LAND TRANSPORT (US DOT):</b>	
DOT Proper Shipping Name:	Not dangerous goods.
DOT Hazard Class:	
UN/NA Number:	
<b>14.1 LAND TRANSPORT (European ADR/RID):</b>	
ADR/RID Shipping Name:	Not dangerous goods.
UN Number:	
Hazard Class:	
<b>14.3 AIR TRANSPORT (ICAO/IATA):</b>	
ICAO/IATA Shipping Name:	Not dangerous goods.
<b>Additional Transport Information:</b>	Transport in accordance with local, state, and federal regulations.



# SAFETY DATA SHEET

## 5-bromo-5-nitro-1,3-Dioxane

Page: 6 of 6

Revision: 12/22/2016

### Section 15. Regulatory Information

#### EPA SARA (Superfund Amendments and Reauthorization Act of 1986) Lists

CAS #	Hazardous Components (Chemical Name)	S. 302 (EHS)	S. 304 RQ	S. 313 (TRI)
30007-47-7	Bromonitrodioxane	No	No	No

CAS #	Hazardous Components (Chemical Name)	Other US EPA or State Lists
30007-47-7	Bromonitrodioxane	CAA HAP,ODC: No; CWA NPDES: No; TSCA: Yes - Inventory; CA PROP.65: No

**Regulatory Information Statement:** This SDS was prepared in accordance with 29 CFR 1910.1200 and Regulation (EC) No.1272/2008.

### Section 16. Other Information

**Revision Date:** 12/22/2016

**Additional Information About This Product:** No data available.

**Company Policy or Disclaimer:** DISCLAIMER: This information is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes.