

CERTIFICATE OF ANALYSIS

Catalog #: 7305 **Lot #** M30892017

Description: HSV-1 Antigen

Viral Strain: F (ATCC)

Source: Grown in Vero Cells
Product manufactured in the United States using only
domestic animal materials processed in U.S. facilities.

Format: Liquid **Exp. Date:** 20 SEP 2022

Protein Concentration: OD 260/280nm Assay:
Acceptable Range: 0.25 – 0.35 mg/mL
Testing Results: 0.33 mg/mL

Activity: ELISA Assay: Acceptable Range: -15% - +50% of Reference Antigen
Testing Results: Pass

Buffer: 0.1 M Glycine Buffer, pH 9.3 – 9.7

Preservative: None

Applications: Each laboratory should determine an optimum working titer for use in its particular application

Storage: Store at $\leq -65^{\circ}\text{C}$. Avoid multiple freeze/thaw cycles

Inactivation: Product inactivated by chemical disruption and verified by Infectivity Testing (Cell Culture).
Acceptable Range: 0 CPE
Testing Results: Pass

Safety Notes: Refer to the appropriate Safety Data Sheet (SDS) for additional information
Since no test method can offer complete assurance that infectious agents are absent, this material should be handled with caution.


Signature

19 October 2017
Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY