



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road

Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: A01278H **Lot #:** 18L33515

Description: hCG Beta \geq 98% Pure
Human Chorionic Gonadotropin beta subunit (hCG-b)

Source: Human Pregnancy Urine

Format: Purified, Lyophilized **Exp. Date:** 18 Nov 2020

Purification: > 99% pure (SDS-PAGE).
< 1% alpha-hCG by ELISA.

Concentration: 4.94 mg/mL (prior to lyophilization) (OD280nm, E^{0.1%} = 0.21).

Buffer: Lyophilized from 30 mM Ammonium Bicarbonate.

Preservative: None

Application: Specific methodologies have not been tested using this product.

Storage: Store lyophilized product at 2-8°C. After reconstitution, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

CORRY THOMAS

26 Jun 2017

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

003