



Meridian
Life Science, Inc.
Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA
Telephone: 901-382-8716
Fax: 901-333-8223
Email: info@meridianlifescience.com
www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: A01451H **Lot #:** 5I26317

Description: Alpha 1 Microglobulin
Alpha 1 Microglobulin, (A1MG) > 98% Pure, Free Form

Source: Human Urine

Format: Purified, Liquid

Purification: > 98% pure (SDS-PAGE).

Concentration: 1.65 mg/mL (A1MG, OD280nm, E = 1.58)

Buffer: 20 mM Tris Buffer, 100 mM Sodium Chloride, pH 8.0 \pm 0.2.

Preservative: 0.05% Sodium Azide

Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store at -20°C. Avoid multiple freeze/thaw cycles.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Brenda Dunn

23 April 2018

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY