



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: A01690H **Lot #:** 1E14617

Description: CA125 Partially Pure Antigen
Cancer Antigen CA125 (Ovarian Cancer) Partially Pure

Source: Human Fluids

Format: Partially Pure, Liquid

Purification:

Purity Ratio:	50.27 kUnits/mL/OD280nm
Known Contaminants:	< 1% Roche Elecsys CA19-9 Assay
	< 1% Roche Elecsys CA15-3 Assay
	< 1% Roche Elecsys 72-4 Assay
	< 1% Roche Elecsys AFP Assay
	< 1% Roche Elecsys CEA Assay
	1% Roche Elecsys Ferritin Assay

Concentration: Specific Protein: 73.90 kUnits/mL Roche (Cobas CA 125 Assay)

Buffer: 50 mM Phosphate, 150 mM Sodium Chloride, pH 7.4 ± 0.2.
Product is 0.2 µm filtered.

Preservative: 0.09% Sodium Azide

Applications: Suitable for use in ELISA and Electrochemiluminescence assay. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store at 2–8°C.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Sincerely, S. Orr

Signature

30 May 2017

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY