



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: A01694D **Lot #:** 6G20617

Description: 25-OH Vitamin D₃ BSA Conjugate
25-Dihydroxy (OH) Vitamin D₃ Bovine Serum Albumin (BSA) Conjugate

Source: Chemically Manufactured

Format: BSA, Lyophilized
Reconstitute 1 mg/vial in 100 µL sterile water.

Purification: ≥ 99% Pure (HPLC)

Concentration: Not Determined

Buffer: Lyophilized from 0.01 M Phosphate Buffered Saline, pH 7.4.

Preservative: None

Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store lyophilized product at 2-8°C. Aliquot and store reconstituted product at -20°C protected from heat and light.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

25 JULY 2017

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY