

5171 Wilfong Road Memphis, TN 38134

Telephone: 901-382-8716 Fax: 901-333-8223

Email: info@meridianlifescience.com www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: B65807B **Lot #:** 4D09115

Description: Goat A' HBsAg (ad/ay)

Goat Antibody to Hepatitis B Surface Antigen (HBsAg), ad/ay Subtype

Biotin Conjugated

Specificity: Reacts with Hepatitis B virus, purified surface antigen.

Host Animal: Goat

Immunogen: Mixture of subtypes ad & ay. Hepatitis B surface antigen purified from human serum.

Format: Biotin, Liquid

Purification: Purified IgG fraction of the antiserum covalently coupled with the N-Hydroxysuccinimide ester of Biotin

under mild conditions to give a high degree of substitution.

Concentration: 4-5 mg/mL (OD280nm, $E^{0.1\%} = 1.4$)

Buffer: 0.01 M Phosphate Buffered Saline, pH 7.2

This product contains no stabilizing proteins.

Preservatives: 0.1% Sodium Azide

Applications: Suitable for use with avidin and streptavidin amplifications systems for ELISA and IEP. Each laboratory

should determine an optimum working titer for use in its particular application. Other applications have not

been tested but use in such assays should not necessarily be excluded.

Storage: Short-term (up to 6 months) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple

freeze/thaw cycles.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

06 January 2017

Date