



Safety Data Sheet  
acc. to OSHA HCS

Printing date 04/25/2017

Reviewed on 12/16/2016

**1 Identification**

- **Product identifier**
- **Trade name:** Endo-Wash Buffer
- **Article number:** D4036-3-6, D4036-3-15, D4036-3-30, D4036-3-60, D4036-3-120, D4036-3-160, D4036-3-240
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Zymo Research Corp.  
17062 Murphy Ave.  
Irvine, CA 92614  
USA  
Phone: 1-949-679-1190 or 1-888-882-9682  
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**  
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

**2 Hazard(s) identification**

- **Classification of the substance or mixture**



GHS02 Flame

Flam. Liq. 2 H225 Highly flammable liquid and vapor.



GHS07

Acute Tox. 4 H302 Harmful if swallowed.  
Skin Irrit. 2 H315 Causes skin irritation.  
Eye Irrit. 2A H319 Causes serious eye irritation.  
STOT SE 3 H336 May cause drowsiness or dizziness.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS02, GHS07
- **Signal word** Danger
- **Hazard-determining components of labeling:**  
guanidinium chloride  
propan-2-ol
- **Hazard statements**  
Highly flammable liquid and vapor.  
Harmful if swallowed.  
Causes skin irritation.  
Causes serious eye irritation.  
May cause drowsiness or dizziness.

(Contd. on page 2)



**Safety Data Sheet**  
*acc. to OSHA HCS*

Printing date 04/25/2017

Reviewed on 12/16/2016

**Trade name: Endo-Wash Buffer**

(Contd. of page 1)

**· Precautionary statements**

- Keep away from heat/sparks/open flames/hot surfaces. No smoking.
- Use explosion-proof electrical/ventilating/lighting/equipment.
- Avoid breathing dust/fume/gas/mist/vapors/spray
- Wear protective gloves / eye protection / face protection.
- Wear protective gloves / eye protection / face protection.
- Ground/bond container and receiving equipment.
- Use only non-sparking tools.
- Take precautionary measures against static discharge.
- Wash thoroughly after handling.
- Do not eat, drink or smoke when using this product.
- Use only outdoors or in a well-ventilated area.
- If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
- If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- Specific treatment (see on this label).
- IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
- IF INHALED: Remove person to fresh air and keep comfortable for breathing.
- If skin irritation occurs: Get medical advice/attention.
- If eye irritation persists: Get medical advice/attention.
- Rinse mouth.
- In case of fire: Use for extinction: CO2, powder or water spray.
- Take off contaminated clothing and wash it before reuse.
- Store in a well-ventilated place. Keep container tightly closed.
- Store in a well-ventilated place. Keep cool.
- Store locked up.
- Dispose of contents/container in accordance with local/regional/national/international regulations.

**· Classification system:**

**· NFPA ratings (scale 0 - 4)**



**· HMIS-ratings (scale 0 - 4)**



**· Other hazards**

**· Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/25/2017

Reviewed on 12/16/2016

Trade name: Endo-Wash Buffer

(Contd. of page 2)

3 Composition/information on ingredients

- Chemical characterization: Mixtures
Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 2 columns: CAS number and chemical name, and a percentage column. Rows include guanidinium chloride (≤50%) and propan-2-ol (≤50%).

4 First-aid measures

- Description of first aid measures
General information: Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.
After inhalation: In case of unconsciousness place patient stably in side position for transportation.
After skin contact: Immediately wash with water and soap and rinse thoroughly.
After eye contact: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
After swallowing: Do not induce vomiting; immediately call for medical help.
Information for doctor:
Most important symptoms and effects, both acute and delayed No further relevant information available.
Indication of any immediate medical attention and special treatment needed No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures Wear protective equipment. Keep unprotected persons away.
Environmental precautions: Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.
Methods and material for containment and cleaning up: Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust). Dispose contaminated material as waste according to item 13. Ensure adequate ventilation.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/25/2017

Reviewed on 12/16/2016

Trade name: Endo-Wash Buffer

(Contd. of page 3)

- Reference to other sections
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include guanidinium chloride (1.4 mg/m3) and propan-2-ol (400 ppm).

PAC-2:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include guanidinium chloride (16 mg/m3) and propan-2-ol (2000\* ppm).

PAC-3:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include guanidinium chloride (94 mg/m3) and propan-2-ol (12000\*\* ppm).

7 Handling and storage

- Handling:
Precautions for safe handling
Information about protection against explosions and fires:
Conditions for safe storage, including any incompatibilities
Storage:
Requirements to be met by storerooms and receptacles:
Information about storage in one common storage facility:
Further information about storage conditions:
Specific end use(s) Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems:
Control parameters
Components with limit values that require monitoring at the workplace:
At this time, the remaining constituent has no known exposure limits.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/25/2017

Reviewed on 12/16/2016

Trade name: Endo-Wash Buffer

(Contd. of page 4)

CAS: 67-63-0 propan-2-ol

Table with 2 columns: PEL, REL, TLV and their corresponding values for propan-2-ol.

Ingredients with biological limit values:

CAS: 67-63-0 propan-2-ol

Table with 2 columns: BEI and its details (Medium, Time, Parameter) for propan-2-ol.

Additional information: The lists that were valid during the creation were used as basis.

Exposure controls

Personal protective equipment:

General protective and hygienic measures:

- Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.

Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/25/2017

Reviewed on 12/16/2016

Trade name: Endo-Wash Buffer

(Contd. of page 5)

Eye protection:



Tightly sealed goggles

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

- Form: Liquid
Color: Clear
Odor: Alcohol-like
Odor threshold: Not determined.

pH-value: Not determined.

Change in condition

- Melting point/Melting range: Undetermined.
Boiling point/Boiling range: 82 °C (180 °F)

Flash point: 13 °C (55 °F)

Flammability (solid, gaseous): Not applicable.

Ignition temperature: 425 °C (797 °F)

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product is not explosive. However, formation of explosive air/vapor mixtures are possible.

Explosion limits:

- Lower: 2.0 Vol %
Upper: 12.0 Vol %

Vapor pressure at 20 °C (68 °F): 43 hPa (32 mm Hg)

- Density: Not determined.
Relative density: Not determined.
Vapor density: Not determined.
Evaporation rate: Not determined.

Solubility in / Miscibility with Water:

Fully miscible.

Partition coefficient (n-octanol/water): Not determined.

Viscosity:

Dynamic: Not determined.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/25/2017

Reviewed on 12/16/2016

Trade name: Endo-Wash Buffer

(Contd. of page 6)

Table with 2 columns: Property (Kinematic, Solvent content, Organic solvents, VOC content, Solids content, Other information) and Value (Not determined, 50.0%, 50.0%, 500.0 g/l / 4.17 lb/gl, 40.0%, No further relevant information available).

10 Stability and reactivity

- Reactivity: No further relevant information available.
Chemical stability
Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
Possibility of hazardous reactions: No dangerous reactions known.
Conditions to avoid: No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:

Table with 3 columns: Route, LD50/LC50, and Value. Includes data for CAS: 50-01-1 guanidinium chloride and CAS: 67-63-0 propan-2-ol.

- Primary irritant effect:
on the skin: Irritant to skin and mucous membranes.
on the eye: Irritating effect.
Sensitization: No sensitizing effects known.
Additional toxicological information:
The product shows the following dangers according to internally approved calculation methods for preparations: Harmful, Irritant.

- Carcinogenic categories

Table with 2 columns: Agency (IARC) and Substance (CAS: 67-63-0 propan-2-ol). Page number 3 is in the bottom right corner.

(Contd. on page 8)



Safety Data Sheet acc. to OSHA HCS

Printing date 04/25/2017

Reviewed on 12/16/2016

Trade name: Endo-Wash Buffer

(Contd. of page 7)

· NTP (National Toxicology Program)

None of the ingredients is listed.

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
· Aquatic toxicity: No further relevant information available.
· Persistence and degradability No further relevant information available.
· Behavior in environmental systems:
· Bioaccumulative potential No further relevant information available.
· Mobility in soil No further relevant information available.
· Additional ecological information:
· General notes:
Water hazard class 1 (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
· Results of PBT and vPvB assessment
· PBT: Not applicable.
· vPvB: Not applicable.
· Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
· Recommendation:
Dispose of contents in accordance with local/regional/national, and international recommendations.
· Uncleaned packagings:
· Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
· Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

- UN-Number
· DOT, IMDG, IATA UN1219
· UN proper shipping name
· DOT Isopropanol mixture
· IMDG, IATA ISOPROPANOL (ISOPROPYL ALCOHOL) mixture

(Contd. on page 9)





Safety Data Sheet
acc. to OSHA HCS

Printing date 04/25/2017

Reviewed on 12/16/2016

Trade name: Endo-Wash Buffer

(Contd. of page 8)

· Transport hazard class(es)

· DOT



· Class 3 Flammable liquids
· Label 3

· IMDG, IATA



· Class 3 Flammable liquids
· Label 3

· Packing group

· DOT, IMDG, IATA II

· Environmental hazards: Not applicable.

· Special precautions for user Warning: Flammable liquids

· Danger code (Kemler): 33

· EMS Number: F-E,S-D

· Stowage Category B

· Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code Not applicable.

· Transport/Additional information:

· DOT

· Quantity limitations On passenger aircraft/rail: 5 L
On cargo aircraft only: 60 L

· IMDG

· Limited quantities (LQ) 1L

· Excepted quantities (EQ) Code: E2
Maximum net quantity per inner packaging: 30 ml
Maximum net quantity per outer packaging: 500 ml

· UN "Model Regulation": UN 1219 ISOPROPANOL MIXTURE, 3, II

US

(Contd. on page 10)



Safety Data Sheet acc. to OSHA HCS

Printing date 04/25/2017

Reviewed on 12/16/2016

Trade name: Endo-Wash Buffer

(Contd. of page 9)

15 Regulatory information

- Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Section 355 (extremely hazardous substances):

None of the ingredients is listed.

Section 313 (Specific toxic chemical listings):

CAS: 67-63-0 propan-2-ol

TSCA (Toxic Substances Control Act):

All ingredients are listed.

Proposition 65

Chemicals known to cause cancer:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

CAS: 67-63-0 propan-2-ol

A4

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms GHS02, GHS07

Signal word Danger

Hazard-determining components of labeling:

guanidinium chloride

propan-2-ol

Hazard statements

Highly flammable liquid and vapor.

Harmful if swallowed.

Causes skin irritation.

Causes serious eye irritation.

May cause drowsiness or dizziness.

Precautionary statements

Keep away from heat/sparks/open flames/hot surfaces. No smoking.

Use explosion-proof electrical/ventilating/lighting/equipment.

(Contd. on page 11)



**Safety Data Sheet**  
*acc. to OSHA HCS*

Printing date 04/25/2017

Reviewed on 12/16/2016

**Trade name: Endo-Wash Buffer**

(Contd. of page 10)

Avoid breathing dust/fume/gas/mist/vapors/spray  
Wear protective gloves / eye protection / face protection.  
Wear protective gloves / eye protection / face protection.  
Ground/bond container and receiving equipment.  
Use only non-sparking tools.  
Take precautionary measures against static discharge.  
Wash thoroughly after handling.  
Do not eat, drink or smoke when using this product.  
Use only outdoors or in a well-ventilated area.  
If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.  
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.  
Continue rinsing.  
Specific treatment (see on this label).  
IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.  
IF INHALED: Remove person to fresh air and keep comfortable for breathing.  
If skin irritation occurs: Get medical advice/attention.  
If eye irritation persists: Get medical advice/attention.  
Rinse mouth.  
In case of fire: Use for extinction: CO2, powder or water spray.  
Take off contaminated clothing and wash it before reuse.  
Store in a well-ventilated place. Keep container tightly closed.  
Store in a well-ventilated place. Keep cool.  
Store locked up.  
Dispose of contents/container in accordance with local/regional/national/international regulations.  
· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

**16 Other information**

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.  
Safety Department  
17062 Murphy Ave.  
Irvine, CA 92614  
USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 04/25/2017 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

(Contd. on page 12)



**ZYMO RESEARCH**

*The Beauty of Science is to Make Things Simple*

Page 12/12

**Safety Data Sheet**  
**acc. to OSHA HCS**

Printing date 04/25/2017

Reviewed on 12/16/2016

**Trade name: Endo-Wash Buffer**

(Contd. of page 11)

NFPA: National Fire Protection Association (USA)  
HMIS: Hazardous Materials Identification System (USA)  
VOC: Volatile Organic Compounds (USA, EU)  
LC50: Lethal concentration, 50 percent  
LD50: Lethal dose, 50 percent  
PBT: Persistent, Bioaccumulative and Toxic  
vPvB: very Persistent and very Bioaccumulative  
NIOSH: National Institute for Occupational Safety  
OSHA: Occupational Safety & Health  
TLV: Threshold Limit Value  
PEL: Permissible Exposure Limit  
REL: Recommended Exposure Limit  
BEI: Biological Exposure Limit  
Flam. Liq. 2: Flammable liquids – Category 2  
Acute Tox. 4: Acute toxicity – Category 4  
Skin Irrit. 2: Skin corrosion/irritation – Category 2  
Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A  
STOT SE 3: Specific target organ toxicity (single exposure) – Category 3

US