

Printing date 02/06/2018

Reviewed on 12/28/2015

Page 1/11

1 Identification	
· Product identifie	
-	quencing Binding Buffer
	D4050-1-14, D4050-1-55, D4050-1-500
· Application of th	ne substance / the mixture Laboratory Reagent
Manufacturer/S Zymo Research (17062 Murphy A Irvine, CA 92614 USA	Corp. ave. 4 9-1190 or 1-888-882-9682
•	artment: Product safety department
• Emergency telep	hone number: usiness hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190
2 Hazard(s) ide	
· Classification of	the substance or mixture
GHS0	5 Corrosion
P	
GHS0 Skin Corr. 1C	5 Corrosion H314 Causes severe skin burns and eye damage. H318 Causes serious eye damage.
GHS0 Skin Corr. 1C Eye Dam. 1	 5 Corrosion H314 Causes severe skin burns and eye damage. H318 Causes serious eye damage. 7 H302 Harmful if swallowed.
GHS0 Skin Corr. 1C Eye Dam. 1 GHS0 Acute Tox. 4 Acute Tox. 4	 5 Corrosion H314 Causes severe skin burns and eye damage. H318 Causes serious eye damage. 7 H302 Harmful if swallowed. H312 Harmful in contact with skin.
GHS0 Skin Corr. 1C Eye Dam. 1 GHS0 Acute Tox. 4 Acute Tox. 4 Acute Tox. 4 Acute Tox. 4	 5 Corrosion H314 Causes severe skin burns and eye damage. H318 Causes serious eye damage. 7 H302 Harmful if swallowed. H312 Harmful in contact with skin. H332 Harmful if inhaled.
GHS0 Skin Corr. 1C Eye Dam. 1 GHS0 Acute Tox. 4 Acute Tox. 4 Acute Tox. 4 Acute Tox. 4 Acute Tox. 4	 5 Corrosion H314 Causes severe skin burns and eye damage. H318 Causes serious eye damage. 7 H302 Harmful if swallowed. H312 Harmful in contact with skin.
GHS0 Skin Corr. 1C Eye Dam. 1 GHS0 Acute Tox. 4 Acute Tox. 4	 5 Corrosion H314 Causes severe skin burns and eye damage. H318 Causes serious eye damage. 7 H302 Harmful if swallowed. H312 Harmful in contact with skin. H332 Harmful if inhaled. 3 H412 Harmful to aquatic life with long lasting effects.
GHS0 Skin Corr. 1C Eye Dam. 1 GHS0 Acute Tox. 4 Acute Tox	 5 Corrosion H314 Causes severe skin burns and eye damage. H318 Causes serious eye damage. 7 H302 Harmful if swallowed. H312 Harmful in contact with skin. H332 Harmful if inhaled. 3 H412 Harmful to aquatic life with long lasting effects. <i>sents</i> The product is classified and labeled according to the Globally Harmonized System (GHS). <i>ms</i> GHS05, GHS07 ger <i>ning components of labeling:</i> cyanate



Printing date 02/06/2018

Trade name: Sequencing Binding Buffer

Reviewed on 12/28/2015

Page 2/11

(Contd. of page 1) · Precautionary statements Do not breathe dusts or mists. Wash thoroughly after handling. Do not eat, drink or smoke when using this product. Use only outdoors or in a well-ventilated area. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. If swallowed: Call a poison center/doctor if you feel unwell. If swallowed: Rinse mouth. Do NOT induce vomiting. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF INHALED: Remove person to fresh air and keep comfortable for breathing. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a poison center/doctor. Specific treatment (see on this label). Take off contaminated clothing and wash it before reuse. Wash contaminated clothing before reuse. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations. · Classification system: · NFPA ratings (scale 0 - 4) Health = 3Fire = 0Reactivity = 0· HMIS-ratings (scale 0 - 4) HEALTH *3 Health = *3FIRE 0 Fire = 0**REACTIVITY** Reactivity = 0· Other hazards · Results of PBT and vPvB assessment · **PBT:** Not applicable. · vPvB: Not applicable. 3 Composition/information on ingredients · Chemical characterization: Mixtures · Description: Mixture of the substances listed below with nonhazardous additions. · Dangerous components: CAS: 593-84-0 guanidinium thiocyanate 70%

(Contd. on page 3)



Printing date 02/06/2018

Reviewed on 12/28/2015

Trade name: Sequencing Binding Buffer

(Contd. of page 2)

Page 3/11

4 First-aid measures

- · Description of first aid measures
- · General information:
- Immediately remove any clothing soiled by the product.

Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

• After inhalation:

Supply fresh air. If required, provide artificial respiration if trained to do so. Keep patient warm. Consult doctor if symptoms persist.

In case of unconsciousness place patient stably in side position for transportation.

- · After skin contact: Immediately wash with water and soap and rinse thoroughly.
- · After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
- After swallowing:
- Rinse mouth
- DO NOT induce vomiting.
- Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

5 Fire-fighting measures

• Extinguishing media

• Suitable extinguishing agents:

CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam. Use fire fighting measures that suit the environment.

· Special hazards arising from the substance or mixture

Products of thermal decomposition of this material would include hydrogen cyanide, ammonia, and oxides of carbon, nitrogen and sulfur.

- Advice for firefighters
- · Protective equipment: Wear self-contained breathing apparatus for fighting fires involving this material

6 Accidental release measures

Personal precautions, protective equipment and emergency procedures
Wear self-contained breathing apparatus for responding to non-incidental release of this material in which there is the potential for inhalation of vapors, mists or sprays
Wear protective equipment. Keep unprotected persons away.
Environmental precautions:
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
Methods and material for containment and cleaning up:
Absorb with liquid-binding material (sand diatomite acid binders universal binders sawdust)

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust). Use neutralizing agent.

Dispose contaminated material as waste according to item 13.

(Contd. on page 4)

US



Printing date 02/06/2018

Reviewed on 12/28/2015

Trade name: Sequencing Binding Buffer (Contd. of page 3) Ensure adequate ventilation. · Reference to other sections See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment. See Section 13 for disposal information. · Protective Action Criteria for Chemicals · PAC-1: CAS: 593-84-0 guanidinium thiocyanate 0.98 mg/m^3 · PAC-2: CAS: 593-84-0 guanidinium thiocyanate 11 mg/m³ · PAC-3: CAS: 593-84-0 guanidinium thiocyanate 65 mg/m³

7 Handling and storage

· Handling:

- · Precautions for safe handling
- Ensure good ventilation/exhaustion at the workplace.
- Prevent formation of aerosols.
- · Information about protection against explosions and fires: No special measures required.
- \cdot Conditions for safe storage, including any incompatibilities
- Store in cool, dry place. Store in well-ventilated location.
- · Storage:
- · Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Do not store together with acids or strong oxidizers
- · Further information about storage conditions: Keep receptacle tightly sealed.
- · *Specific end use(s)* Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
- · Control parameters

Work under a chemical fume hood when using this product. Ensure eyewash station and safety showers are readily accessible.

· Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls

The appropriate protective equipment under anticipated circumstances of use include lab-coat, safety glasses with side-shields and gloves.

(Contd. on page 5)

Page 4/11

US



Printing date 02/06/2018

Reviewed on 12/28/2015

Trade name: Sequencing Binding Buffer

· Personal protective equipment:

· General protective and hygienic measures:

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

· Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· Protection of hands:



The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

• Eye protection:



Tightly sealed goggles

9 Physical and chemical properties

· Information on basic physical and · General Information	chemical properties	
· Appearance:		
Form:	Liquid	
Color:	Clear	
· Odor:	Mild	
· Odor threshold:	Not determined.	
· pH-value:	Not determined.	
· Change in condition		
Melting point/Melting range:	Undetermined.	
Boiling point/Boiling range:	Undetermined.	
		(Contd. on page 6

(Contd. on page 6)

(Contd. of page 4)



Printing date 02/06/2018

Reviewed on 12/28/2015

Trade name: Sequencing Binding Buffer

	(Contd. of p	age 5
· Flash point:	Not applicable.	
· Flammability (solid, gaseous):	Not applicable.	
· Ignition temperature:		
Decomposition temperature:	Not determined.	
· Auto igniting:	Product is not selfigniting.	
• Danger of explosion:	Product does not present an explosion hazard.	
· Explosion limits:		
Lower:	Not Applicable	
Upper:	Not Applicable	
· Vapor pressure:	Not determined.	
· Density:	Not determined.	
· Relative density	Not determined.	
· Vapor density	Not determined.	
· Evaporation rate	Not determined.	
· Solubility in / Miscibility with		
Water:	Fully miscible.	
· Partition coefficient (n-octanol/wat	ter): Not determined.	
· Viscosity:		
Dynamic:	Not determined.	
Kinematic:	Not determined.	
· Solvent content:		
VOC content:	0.00 %	
	0.0 g/l / 0.00 lb/gl	
Solids content:	70.0 %	
• Other information	No further relevant information available.	

10 Stability and reactivity

- · Reactivity No further relevant information available.
- · Chemical stability This product is normally stable under anticipated circumstances of use and storage.
- *Thermal decomposition / conditions to be avoided:* Products of thermal decomposition of this material would include hydrogen cyanide, ammonia, and oxides of carbon nitrogen and sulfur.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid Avoid exposing product to extreme temperatures or incompatible chemicals
- · Incompatible materials: Acids and strong oxidizers

(Contd. on page 7)

US



Printing date 02/06/2018

Reviewed on 12/28/2015

Trade name: Sequencing Binding Buffer

· Hazardous decomposition products:

Product will not undergo self-decomposition, so no such products will be generated.

11 Toxicological information

- · Information on toxicological effects
- Acute toxicity:
- · LD/LC50 values that are relevant for classification:

CAS: 593-84-0 guanidinium thiocyanate

Oral LD50 593 mg/kg (rat)

- · Primary irritant effect:
- · on the skin: Caustic effect on skin and mucous membranes.
- · on the eye: Strong caustic effect.
- · Sensitization: No sensitizing effects known.
- \cdot Additional toxicological information:

The product shows the following dangers according to internally approved calculation methods for preparations: Harmful

Corrosive

Swallowing will lead to a strong caustic effect on mouth and throat and to the danger of perforation of esophagus and stomach.

· Carcinogenic categories

· IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

· Toxicity

· Aquatic toxicity:

CAS: 593-84-0 guanidinium thiocyanate

EC50 42.4 mg/kg (daphnia)

· Persistence and degradability No further relevant information available.

· Behavior in environmental systems:

- · *Bioaccumulative potential* No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:

· General notes:

Water hazard class 2 (Self-assessment): hazardous for water

Do not allow product to reach ground water, water course or sewage system.

(Contd. on page 8)

Page 7/11

(Contd. of page 6)



Printing date 02/06/2018

Reviewed on 12/28/2015

Trade name: Sequencing Binding Buffer

Must not reach bodies of water or drainage ditch undiluted or unneutralized. Danger to drinking water if even small quantities leak into the ground.

· Results of PBT and vPvB assessment

- *PBT*: Not applicable.
- *vPvB*: Not applicable.
- · Other adverse effects No further relevant information available.

13 Disposal considerations

· Waste treatment methods

· Recommendation:

Dispose of contents in accordance with local/regional/national, and international recommendations.

- · Uncleaned packagings:
- · Recommendation:

Dispose of container in acoordance with local/regional/national and international recommendations.

• *Recommended cleansing agent:* Water, if necessary with cleansing agents.

UN1760
Corrosive liquids, n.o.s. (guanidinium thiocyanate)
CORROSIVE LIQUID, N.O.S. (guanidinium thiocyanate)
8 Corrosive substances
8
8 Corrosive substances
8

Page 8/11

(Contd. of page 7)



Printing date 02/06/2018

Reviewed on 12/28/2015

Trade name: Sequencing Binding Buffer

	(Contd. of page
· Environmental hazards:	Not applicable.
· Special precautions for user	Warning: Corrosive substances
· Danger code (Kemler):	80
· EMS Number:	F-A,S-B
· Stowage Category	А
· Stowage Code	SW2 Clear of living quarters.
• Transport in bulk according to Annex	II of
MARPOL73/78 and the IBC Code	Not applicable.
· Transport/Additional information:	
·DOT	
• Quantity limitations	On passenger aircraft/rail: 5 L
~ ·	On cargo aircraft only: 60 L
· IMDG	
· Limited quantities (LQ)	5L
· Excepted quantities $(\widetilde{E}Q)$	Code: E1
	Maximum net quantity per inner packaging: 30 ml
	Maximum net quantity per outer packaging: 1000 ml
· UN "Model Regulation":	UN 1760 CORROSIVE LIQUIDS, N.O.S. (GUANIDINIUM THIOCYANATE), 8, III

15 Regulatory information

 \cdot Safety, health and environmental regulations/legislation specific for the substance or mixture \cdot Sara

• Section 355 (extremely hazardous substances):
None of the ingredients is listed.
· Section 313 (Specific toxic chemical listings):
None of the ingredients is listed.
· TSCA (Toxic Substances Control Act):
All ingredients are listed.
· Proposition 65
· Chemicals known to cause cancer:
None of the ingredients is listed.
· Chemicals known to cause reproductive toxicity for females:
None of the ingredients is listed.
· Chemicals known to cause reproductive toxicity for males:
None of the ingredients is listed.

(Contd. on page 10)

Page 9/11

US



Page 10/11

(Contd. of page 9)

Safety Data Sheet acc. to OSHA HCS

Printing date 02/06/2018

Reviewed on 12/28/2015

Trade name: Sequencing Binding Buffer

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

• TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

- · GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).
- · Hazard pictograms GHS05, GHS07
- · Signal word Danger
- · Hazard-determining components of labeling: guanidinium thiocyanate · Hazard statements Harmful if swallowed, in contact with skin or if inhaled. Causes severe skin burns and eye damage. Harmful to aquatic life with long lasting effects. · Precautionary statements Do not breathe dusts or mists. Wash thoroughly after handling. Do not eat, drink or smoke when using this product. Use only outdoors or in a well-ventilated area. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. If swallowed: Call a poison center/doctor if you feel unwell. If swallowed: Rinse mouth. Do NOT induce vomiting. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF INHALED: Remove person to fresh air and keep comfortable for breathing. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a poison center/doctor. Specific treatment (see on this label). Take off contaminated clothing and wash it before reuse. Wash contaminated clothing before reuse. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations. · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

(Contd. on page 11)

____US



Printing date 02/06/2018

Reviewed on 12/28/2015

Trade name: Sequencing Binding Buffer

(Contd. of page 10)
· Department issuing SDS:
Zymo Research Corp.
Safety Department
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
· <i>Contact:</i> sds@zymoresearch.com
· Date of preparation / last revision 02/06/2018 / -
· Abbreviations and acronyms:
ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage
of Dangerous Goods by Road)
IMDG: International Maritime Code for Dangerous Goods
DOT: US Department of Transportation
IATA: International Air Transport Association
ACGIH: American Conference of Governmental Industrial Hygienists
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA) HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
LC50: Lethal concentration, 50 percent
LD50: Lethal dose, 50 percent
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety
OSHA: Occupational Safety & Health
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit
Acute Tox. 4: Acute toxicity – Category 4
Skin Corr. 1C: Skin corrosion/irritation – Category 1C
Eye Dam. 1: Serious eye damage/eye irritation – Category 1 Aquatic Chronic 3: Hazardous to the aquatic environment - long-term aquatic hazard – Category 3
Aquate Chrome 5. Hazardous to the aquate environment - tong-term aquate hazard – Category 5

Page 11/11