



**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 10/11/2017

Reviewed on 02/02/2016

**1 Identification**

- **Product identifier**
- **Trade name:** ZymoPURE Wash 2 (Concentrate)
- **Article number:** D4200-6-10, D4200-6-12, D4200-6-23, D4200-6-28
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Zymo Research Corp.  
17062 Murphy Ave.  
Irvine, CA 92614  
USA  
Phone: 1-949-679-1190 or 1-888-882-9682  
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**  
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

**2 Hazard(s) identification**

- **Classification of the substance or mixture**



GHS02 Flame

Flam. Liq. 2 H225 Highly flammable liquid and vapor.



GHS07

Eye Irrit. 2A H319 Causes serious eye irritation.

STOT SE 3 H336 May cause drowsiness or dizziness.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS02, GHS07
- **Signal word** Danger
- **Hazard-determining components of labeling:**  
propan-2-ol
- **Hazard statements**  
Highly flammable liquid and vapor.  
Causes serious eye irritation.  
May cause drowsiness or dizziness.
- **Precautionary statements**  
Keep away from heat/sparks/open flames/hot surfaces. No smoking.  
Use explosion-proof electrical/ventilating/lighting/equipment.  
Avoid breathing dust/fume/gas/mist/vapors/spray  
Wear protective gloves / eye protection / face protection.  
Ground/bond container and receiving equipment.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 10/11/2017

Reviewed on 02/02/2016

Trade name: ZymoPURE Wash 2 (Concentrate)

(Contd. of page 1)

Use only non-sparking tools.
Take precautionary measures against static discharge.
Wash thoroughly after handling.
Use only outdoors or in a well-ventilated area.
If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
Continue rinsing.
IF INHALED: Remove person to fresh air and keep comfortable for breathing.
Call a POISON CENTER/doctor if you feel unwell.
If eye irritation persists: Get medical advice/attention.
In case of fire: Use for extinction: CO2, powder or water spray.
Store in a well-ventilated place. Keep container tightly closed.
Store in a well-ventilated place. Keep cool.
Store locked up.
Dispose of contents/container in accordance with local/regional/national/international regulations.

Classification system:

NFPA ratings (scale 0 - 4)



HMIS-ratings (scale 0 - 4)



Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

3 Composition/information on ingredients

Chemical characterization: Mixtures

Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 2 columns: CAS: 67-63-0, propan-2-ol, and ≤100%

4 First-aid measures

Description of first aid measures

After inhalation: Supply fresh air; consult doctor in case of complaints.

After skin contact: Generally the product does not irritate the skin.

After eye contact:

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 10/11/2017

Reviewed on 02/02/2016

Trade name: ZymoPURE Wash 2 (Concentrate)

(Contd. of page 2)

- After swallowing: Do not induce vomiting; immediately call for medical help.
Information for doctor:
Most important symptoms and effects, both acute and delayed No further relevant information available.
Indication of any immediate medical attention and special treatment needed
No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures
Wear protective equipment. Keep unprotected persons away.
Environmental precautions:
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
Methods and material for containment and cleaning up:
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Ensure adequate ventilation.
Reference to other sections
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
Protective Action Criteria for Chemicals

PAC-1:

CAS: 67-63-0 propan-2-ol 400 ppm

PAC-2:

CAS: 67-63-0 propan-2-ol 2000\* ppm

PAC-3:

CAS: 67-63-0 propan-2-ol 12000\*\* ppm

7 Handling and storage

- Handling:
Precautions for safe handling
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.
Information about protection against explosions and fires:
Keep ignition sources away - Do not smoke.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 10/11/2017

Reviewed on 02/02/2016

Trade name: ZymoPURE Wash 2 (Concentrate)

(Contd. of page 3)

Protect against electrostatic charges.

Conditions for safe storage, including any incompatibilities

Storage:

Requirements to be met by storerooms and receptacles: Store in a cool location.

Information about storage in one common storage facility: Not required.

Further information about storage conditions:

Keep receptacle tightly sealed.

Store in cool, dry conditions in well sealed receptacles.

Specific end use(s) No further relevant information available.

8 Exposure controls/personal protection

Additional information about design of technical systems: No further data; see item 7.

Control parameters

Components with limit values that require monitoring at the workplace:

CAS: 67-63-0 propan-2-ol

PEL Long-term value: 980 mg/m³, 400 ppm

REL Short-term value: 1225 mg/m³, 500 ppm

Long-term value: 980 mg/m³, 400 ppm

TLV Short-term value: 984 mg/m³, 400 ppm

Long-term value: 492 mg/m³, 200 ppm

BEI

Ingredients with biological limit values:

CAS: 67-63-0 propan-2-ol

BEI 40 mg/L

Medium: urine

Time: end of shift at end of workweek

Parameter: Acetone (background, nonspecific)

Additional information: The lists that were valid during the creation were used as basis.

Exposure controls

Personal protective equipment:

General protective and hygienic measures:

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes.

Avoid contact with the eyes and skin.

Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

(Contd. on page 5)



**Safety Data Sheet**  
*acc. to OSHA HCS*

Printing date 10/11/2017

Reviewed on 02/02/2016

**Trade name: ZymoPURE Wash 2 (Concentrate)**

(Contd. of page 4)

**· Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

**· Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

**· Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

**· Eye protection:**



Tightly sealed goggles

**9 Physical and chemical properties**

**· Information on basic physical and chemical properties**

**· General Information**

**· Appearance:**

<b>Form:</b>	Liquid
<b>Color:</b>	Clear
<b>· Odor:</b>	Alcohol-like
<b>· Odor threshold:</b>	Not determined.

**· pH-value:** Not determined.

**· Change in condition**

<b>Melting point/Melting range:</b>	Undetermined.
<b>Boiling point/Boiling range:</b>	82 °C (180 °F)

**· Flash point:** 13 °C (55 °F)

**· Flammability (solid, gaseous):** Not applicable.

**· Ignition temperature:** 425 °C (797 °F)

**· Decomposition temperature:** Not determined.

**· Auto igniting:** Product is not selfigniting.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 10/11/2017

Reviewed on 02/02/2016

Trade name: ZymoPURE Wash 2 (Concentrate)

(Contd. of page 5)

Table with hazard information including: Danger of explosion, Explosion limits (Lower: 2.0 Vol %, Upper: 12.0 Vol %), Vapor pressure at 20 °C (68 °F): 43 hPa (32 mm Hg), Density, Relative density, Vapor density, Evaporation rate, Solubility in / Miscibility with Water: Fully miscible, Partition coefficient (n-octanol/water): Not determined, Viscosity (Dynamic, Kinematic), Solvent content (Organic solvents: 80.0 %, VOC content: 80.0 %), Other information: No further relevant information available.

10 Stability and reactivity

- Reactivity: No further relevant information available.
Chemical stability
Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
Possibility of hazardous reactions: No dangerous reactions known.
Conditions to avoid: No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:

LD/LC50 values that are relevant for classification:

CAS: 67-63-0 propan-2-ol

Table with 3 columns: Route, LD50, and Toxicity Data. Rows: Oral (LD50 5045 mg/kg (rat)), Dermal (LD50 12800 mg/kg (rabbit))

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 10/11/2017

Reviewed on 02/02/2016

Trade name: ZymoPURE Wash 2 (Concentrate)

(Contd. of page 6)

Inhalative LC50/4 h 30 mg/l (rat)

Primary irritant effect:

- on the skin: No irritant effect.
on the eye: Irritating effect.
Sensitization: No sensitizing effects known.
Additional toxicological information:

The product shows the following dangers according to internally approved calculation methods for preparations: Irritant

Carcinogenic categories

IARC (International Agency for Research on Cancer)

CAS: 67-63-0 propan-2-ol

3

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes:
Water hazard class 1 (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation:
Dispose of contents in accordance with local/regional/national, and international recommendations.
Uncleaned packagings:
Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.

(Contd. on page 8)



Safety Data Sheet
acc. to OSHA HCS

Printing date 10/11/2017

Reviewed on 02/02/2016

Trade name: ZymoPURE Wash 2 (Concentrate)

(Contd. of page 7)

Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Property and Value. Rows include UN-Number (UN1219), UN proper shipping name (Isopropanol solution), Transport hazard class(es) (3 Flammable liquids), IMDG, IATA (3 Flammable liquids), Packing group (II), Environmental hazards (Not applicable), Special precautions for user (Warning: Flammable liquids), and Transport/Additional information (Quantity limitations).

(Contd. on page 9)





Safety Data Sheet
acc. to OSHA HCS

Printing date 10/11/2017

Reviewed on 02/02/2016

Trade name: ZymoPURE Wash 2 (Concentrate)

(Contd. of page 8)

- IMDG
Limited quantities (LQ) 1L
Excepted quantities (EQ) Code: E2
Maximum net quantity per inner packaging: 30 ml
Maximum net quantity per outer packaging: 500 ml
UN "Model Regulation": UN 1219 ISOPROPANOL SOLUTION, 3, II

15 Regulatory information

- Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Section 355 (extremely hazardous substances):

None of the ingredients is listed.

Section 313 (Specific toxic chemical listings):

All ingredients are listed.

TSCA (Toxic Substances Control Act):

All ingredients are listed.

Proposition 65

Chemicals known to cause cancer:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

CAS: 67-63-0 propan-2-ol

A4

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms GHS02, GHS07

Signal word Danger

Hazard-determining components of labeling:

propan-2-ol

(Contd. on page 10)

**Safety Data Sheet**  
*acc. to OSHA HCS*

Printing date 10/11/2017

Reviewed on 02/02/2016

**Trade name: ZymoPURE Wash 2 (Concentrate)**

(Contd. of page 9)

**· Hazard statements**

Highly flammable liquid and vapor.  
Causes serious eye irritation.  
May cause drowsiness or dizziness.

**· Precautionary statements**

Keep away from heat/sparks/open flames/hot surfaces. No smoking.  
Use explosion-proof electrical/ventilating/lighting/equipment.  
Avoid breathing dust/fume/gas/mist/vapors/spray  
Wear protective gloves / eye protection / face protection.  
Ground/bond container and receiving equipment.  
Use only non-sparking tools.  
Take precautionary measures against static discharge.  
Wash thoroughly after handling.  
Use only outdoors or in a well-ventilated area.  
If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.  
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
IF INHALED: Remove person to fresh air and keep comfortable for breathing.  
Call a POISON CENTER/doctor if you feel unwell.  
If eye irritation persists: Get medical advice/attention.  
In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.  
Store in a well-ventilated place. Keep container tightly closed.  
Store in a well-ventilated place. Keep cool.  
Store locked up.  
Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

**16 Other information**

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

**· Department issuing SDS:**

Zymo Research Corp.  
Safety Department  
17062 Murphy Ave.  
Irvine, CA 92614  
USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 10/11/2017 / -

**· Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
DOT: US Department of Transportation  
IATA: International Air Transport Association  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service (division of the American Chemical Society)

(Contd. on page 11)



**ZYMO RESEARCH**

*The Beauty of Science is to Make Things Simple*

Page 11/11

**Safety Data Sheet**  
**acc. to OSHA HCS**

Printing date 10/11/2017

Reviewed on 02/02/2016

**Trade name: ZymoPURE Wash 2 (Concentrate)**

(Contd. of page 10)

LC50: Lethal concentration, 50 percent  
LD50: Lethal dose, 50 percent  
PBT: Persistent, Bioaccumulative and Toxic  
vPvB: very Persistent and very Bioaccumulative  
TLV: Threshold Limit Value  
PEL: Permissible Exposure Limit  
REL: Recommended Exposure Limit  
Flam. Liq. 2: Flammable liquids – Category 2  
Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A  
STOT SE 3: Specific target organ toxicity (single exposure) – Category 3

US