

Printing date 12/21/2018

Reviewed on 03/22/2017

Page 1/11

Printing date 12/21/2018	<i>Kevlewea on 05/22/2017</i>
1 Identification	
· Product identifier	
• Trade name: DNA Prep Buffer	
 Article number: D7001-2-10, D7001-2-25, D7001-2-200 Application of the substance / the mixture Laboratory Reagent 	
 Details of the supplier of the safety data sheet Manufacturer/Supplier: Zymo Research Corp. 17062 Murphy Ave. Irvine, CA 92614 USA Phone: 1-949-679-1190 or 1-888-882-9682 sds@zymoresearch.com 	
· Information department: Product safety department	
• <i>Emergency telephone number:</i> During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (9-	49) 679 1190
2 Hazard(s) identification	
· Classification of the substance or mixture	
GHS05 Corrosion	
Skin Corr. 1C H314 Causes severe skin burns and eye damage.	
Eye Dam. 1 H318 Causes serious eye damage.	
Acute Tox. 4 H302 Harmful if swallowed.	
Aquatic Chronic 3 H412 Harmful to aquatic life with long lasting effects.	
 Label elements GHS label elements The product is classified and labeled according to the Hazard pictograms GHS05, GHS07 Signal word Danger 	Globally Harmonized System (GHS).
• <i>Hazard-determining components of labeling:</i> guanidinium thiocyanate • <i>Hazard statements</i>	
Harmful if swallowed.	
Causes severe skin burns and eye damage. Harmful to aquatic life with long lasting effects.	
· Precautionary statements	
If medical advice is needed, have product container or label at hand.	
Keep out of reach of children. Read label before use.	
	(Contd. on page 2)



Printing date 12/21/2018

Reviewed on 03/22/2017

Trade name: DNA Prep Buffer (Contd. of page 1) Do not breathe dusts or mists. Wash thoroughly after handling. Do not eat, drink or smoke when using this product. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. If swallowed: Call a poison center/doctor if you feel unwell. If swallowed: Rinse mouth. Do NOT induce vomiting. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF INHALED: Remove person to fresh air and keep comfortable for breathing. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a poison center/doctor. Specific treatment (see on this label). Wash contaminated clothing before reuse. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations. · Classification system: · NFPA ratings (scale 0 - 4) Health = 3Fire = 1Reactivity = 0· HMIS-ratings (scale 0 - 4) HEALTH Health = *3*3 FIRE Fire = 11 Reactivity = 0REACTIVITY 0 · Other hazards · Results of PBT and vPvB assessment · **PBT**: Not applicable. · vPvB: Not applicable. **3** Composition/information on ingredients · Chemical characterization: Mixtures · Description: Mixture of the substances listed below with nonhazardous additions.

 Dangerous compon 	ients:	
CAS: 593-84-0 gua	anidinium thiocyanate	50%
CAS: 56-81-5 gly	rcerol	50%

(Contd. on page 3)

Page 2/11



Printing date 12/21/2018

Trade name: DNA Prep Buffer

4 First-aid measures

- · Description of first aid measures
- · General information:
- Immediately remove any clothing soiled by the product.

Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

• After inhalation:

Supply fresh air. If required, provide artificial respiration if trained to do so. Keep patient warm. Consult doctor if symptoms persist.

In case of unconsciousness place patient stably in side position for transportation.

· After skin contact: Immediately wash with water and soap and rinse thoroughly.

• After eye contact:

Wash eyes immediately, for at least 15 minutes, with large amounts of water, holding upper and lower lids open. Remove contact lenses, if present and it is easy to do so. Get medical attention immediately.

• After swallowing:

Rinse mouth

DO NOT induce vomiting.

• Information for doctor:

• Most important symptoms and effects, both acute and delayed No further relevant information available. • Indication of any immediate medical attention and special treatment needed

No further relevant information available.

5 Fire-fighting measures

- · Extinguishing media
- · Suitable extinguishing agents: Use fire fighting measures that suit the environment.
- \cdot Special hazards arising from the substance or mixture

Products of thermal decomposition of this material would include hydrogen cyanide, ammonia, and oxides of carbon, nitrogen and sulfur.

- Advice for firefighters
- · Protective equipment: Wear self-contained breathing apparatus for fighting fires involving this material

6 Accidental release measures

*Personal precautions, protective equipment and emergency procedures*Wear self-contained breathing apparatus for responding to non-incidental release of this material in which there is the potential for inhalation of vapors, mists or sprays *Wear* protective equipment. Keep unprotected persons away. *Environmental precautions:*Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water. *Methods and material for containment and cleaning up:*Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Use neutralizing agent. Dispose contaminated material as waste according to item 13.

(Contd. on page 4)

Reviewed on 03/22/2017

(Contd. of page 2)

Page 3/11



Printing date 12/21/2018

Reviewed on 03/22/2017

Page 4/11

Trade name: DNA Prep Buffer		
		(Contd. of page 3)
Ensure adequate		
· Reference to other		
	or information on safe handling.	
	or information on personal protection equipment.	
	for disposal information.	
· Protective Actio	on Criteria for Chemicals	
· PAC-1:		
CAS: 593-84-0	guanidinium thiocyanate	0.98 mg/m ³
CAS: 56-81-5	glycerol	45 mg/m ³
· PAC-2:		
CAS: 593-84-0	guanidinium thiocyanate	11 mg/m ³
CAS: 56-81-5	glycerol	180 mg/m ³
· PAC-3:		
CAS: 593-84-0	guanidinium thiocyanate	65 mg/m ³
CAS: 56-81-5	glycerol	1,100 mg/m ³

7 Handling and storage

- · Handling:
- *Precautions for safe handling* Ensure good ventilation/exhaustion at the workplace. Prevent formation of aerosols.
- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- Store in cool, dry place. Store in well-ventilated location.
- · Storage:
- Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Do not store together with acids or strong oxidizers
- · Further information about storage conditions: Keep receptacle tightly sealed.
- · Specific end use(s) Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
- · Control parameters

Work under a chemical fume hood when using this product. Ensure eyewash station and safety showers are readily accessible.

· Components with limit values that require monitoring at the workplace:

The following constituent is the only constituent of the product which has a PEL, TLV or other recommended exposure limit.

At this time, the remaining constituent has no known exposure limits.

(Contd. on page 5)

US



Printing date 12/21/2018

Trade name: DNA Prep Buffer

Reviewed on 03/22/2017

(Contd. of page 4)

CAS: 56-81-5 glycerol PEL Long-term value: 15* 5** mg/m³

mist; *total dust **respirable fraction

TLV TLV withdrawn-insufficient data human occup. exp.

• Additional information: The lists that were valid during the creation were used as basis.

· Exposure controls

The appropriate protective equipment under anticipated circumstances of use include lab-coat, safety glasses with side-shields and gloves.

· Personal protective equipment:

· General protective and hygienic measures:

Keep away from foodstuffs, beverages and feed. Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

· Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.


Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation • *Material of gloves*

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

• Eye protection:



(Contd. on page 6)



Printing date 12/21/2018

Trade name: DNA Prep Buffer

9 Physical and chemical properties

Information on basic physical and c	chemical properties
· General Information	
· Appearance: Form:	I :: d
	Liquid Colorless
Color: • Odor:	No information available
· Odor: · Odor threshold:	Not determined.
· pH-value:	Not determined.
-	
· Change in condition	TT-1 (complex)
Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	290 °C (554 °F)
· Flash point:	160 °C (320 °F)
· Flammability (solid, gaseous):	Not applicable.
· Ignition temperature:	400 °C (752 °F)
• Decomposition temperature:	Not determined.
· Auto igniting:	Product is not selfigniting.
Danger of explosion:	Product does not present an explosion hazard.
· Explosion limits:	
Lower:	0.9 Vol %
Upper:	0.0 Vol %
· Vapor pressure at 20 °C (68 °F):	0.1 hPa (0.1 mm Hg)
• Density at 20 •C (68 •F):	644.63 g/cm ³ (5,379.43735 lbs/gal)
· Relative density	Not determined.
· Vapor density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with	
Water:	Fully miscible.
· Partition coefficient (n-octanol/wate	er): Not determined.
· Viscosity:	
Dynamic:	Not determined.
Kinematic:	Not determined.
· Solvent content:	
Organic solvents:	50.0 %
VOC content:	0.00 %
	0.0 g/l / 0.00 lb/gal
Solids content:	50.0 %
	(Contd. on pag

Page 6/11

(Contd. of page 5)

US

Reviewed on 03/22/2017



Printing date 12/21/2018

Reviewed on 03/22/2017

Trade name: DNA Prep Buffer

• Other information

No further relevant information available.

10 Stability and reactivity

- *Reactivity* No further relevant information available.
- · Chemical stability This product is normally stable under anticipated circumstances of use and storage.
- *Thermal decomposition / conditions to be avoided:* Products of thermal decomposition of this material would include hydrogen cyanide, ammonia, and oxides of carbon nitrogen and sulfur.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid Avoid exposing product to extreme temperatures or incompatible chemicals
- · Incompatible materials: Acids and strong oxidizers
- *Hazardous decomposition products:* Product will not undergo self-decomposition, so no such products will be generated.

11 Toxicological information

- · Information on toxicological effects
- Acute toxicity:
- · LD/LC50 values that are relevant for classification:

CAS: 593-84-0 guanidinium thiocyanate

Oral LD50 593 mg/kg (rat)

· Primary irritant effect:

- · on the skin: Caustic effect on skin and mucous membranes.
- · on the eye: Strong caustic effect.
- · Sensitization: No sensitizing effects known.
- · Additional toxicological information:

The product shows the following dangers according to internally approved calculation methods for preparations: Harmful

Corrosive

Swallowing will lead to a strong caustic effect on mouth and throat and to the danger of perforation of esophagus and stomach.

· Carcinogenic categories

· IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

(Contd. on page 8)

Page 7/11

(Contd. of page 6)



Printing date 12/21/2018

Reviewed on 03/22/2017

Trade name: DNA Prep Buffer

(Contd. of page 7)

12 Ecological information

· Toxicity

· Aquatic toxicity:

CAS: 593-84-0 guanidinium thiocyanate

- EC50 42.4 mg/kg (daphnia)
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · *Bioaccumulative potential* No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes:

Water hazard class 2 (Self-assessment): hazardous for water

Do not allow product to reach ground water, water course or sewage system. Must not reach bodies of water or drainage ditch undiluted or unneutralized.

Danger to drinking water if even small quantities leak into the ground.

- · Results of PBT and vPvB assessment
- · *PBT*: Not applicable.
- · *vPvB*: Not applicable.
- · Other adverse effects No further relevant information available.

13 Disposal considerations

- · Waste treatment methods
- · Recommendation:

Dispose of contents in accordance with local/regional/national, and international recommendations.

- · Uncleaned packagings:
- · Recommendation:
- Dispose of container in acoordance with local/regional/national and international recommendations.
- · Recommended cleansing agent: Water, if necessary with cleansing agents.

UN-Number	
DOT, IMDG, IATA	UN1760
UN proper shipping name	
DOT	Corrosive liquids, n.o.s. (guanidinium thiocyanate)
IMDG, IATA	CORROSIVE LIQUID, N.O.S. (guanidinium thiocyanate)



Printing date 12/21/2018

Reviewed on 03/22/2017

	(Contd. of pag
Transport hazard class(es)	
DOT	
• Class	8 Corrosive substances
Label	8
IMDG, IATA	
Class	8 Corrosive substances
Label	8
Packing group DOT, IMDG, IATA	III
Environmental hazards:	Not applicable.
Special precautions for user	Warning: Corrosive substances
Danger code (Kemler):	80
EMS Number:	F-A,S-B
Stowage Category	А
Stowage Code	SW2 Clear of living quarters.
Transport in bulk according to Annex A MARPOL73/78 and the IBC Code	II of Not applicable.
Transport/Additional information:	
DOT	
Quantity limitations	On passenger aircraft/rail: 5 L On cargo aircraft only: 60 L
IMDG	
Limited quantities (LQ)	5L
Excepted quantities (\widetilde{EQ})	Code: E1
	Maximum net quantity per inner packaging: 30 ml Maximum net quantity per outer packaging: 1000 ml
UN "Model Regulation":	UN 1760 CORROSIVE LIQUIDS, N.O.S. (GUANIDINIUM THIOCYANATE), 8, III

(Contd. on page 10)

Page 9/11



Printing date 12/21/2018

Trade name: DNA Prep Buffer

Reviewed on 03/22/2017

(Contd. of page 9)

Page 10/11

15 Regulatory information

 \cdot Safety, health and environmental regulations/legislation specific for the substance or mixture \cdot Sara

· Section 355 (extremely hazardous substances):

None of the ingredients is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

· TSCA (Toxic Substances Control Act):

All ingredients are listed.

· Proposition 65

· Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

• *GHS label elements* The product is classified and labeled according to the Globally Harmonized System (GHS). • *Hazard pictograms* GHS05, GHS07

· Signal word Danger

- Hazard-determining components of labeling: guanidinium thiocyanate
 Hazard statements
 Harmful if swallowed.
 Causes severe skin burns and eye damage.
 Harmful to aquatic life with long lasting effects.
 Precautionary statements
 If medical advice is needed, have product container or label at hand. Keep out of reach of children.
- Reep out of reach of chil

Read label before use.

Do not breathe dusts or mists. Wash thoroughly after handling.



Printing date 12/21/2018

Reviewed on 03/22/2017

Trade name: DNA Prep Buffer

(Contd. of page 10) Do not eat, drink or smoke when using this product. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. If swallowed: Call a poison center/doctor if you feel unwell. If swallowed: Rinse mouth. Do NOT induce vomiting. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF INHALED: Remove person to fresh air and keep comfortable for breathing. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a poison center/doctor. Specific treatment (see on this label). Wash contaminated clothing before reuse. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations. · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS: Zymo Research Corp. Safety Department 17062 Murphy Ave. Irvine, CA 92614 USA Phone: 1-949-679-1190 or 1-888-882-9682 · Contact: sds@zymoresearch.com · Date of preparation / last revision 12/21/2018 / -· Abbreviations and acronyms: IMDG: International Maritime Code for Dangerous Goods DOT: US Department of Transportation IATA: International Air Transport Association ACGIH: American Conference of Governmental Industrial Hygienists EINECS: European Inventory of Existing Commercial Chemical Substances ELINCS: European List of Notified Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society) NFPA: National Fire Protection Association (USA) HMIS: Hazardous Materials Identification System (USA) VOC: Volatile Organic Compounds (USA, EU) LC50: Lethal concentration, 50 percent LD50: Lethal dose, 50 percent PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative NIOSH: National Institute for Occupational Safety OSHA: Occupational Safety & Health TLV: Threshold Limit Value PEL: Permissible Exposure Limit **REL: Recommended Exposure Limit** Acute Tox. 4: Acute toxicity - Category 4 Skin Corr. 1C: Skin corrosion/irritation - Category 1C Eye Dam. 1: Serious eye damage/eye irritation - Category 1 Aquatic Chronic 3: Hazardous to the aquatic environment - long-term aquatic hazard - Category 3

Page 11/11