

Printing date 12/12/2016

Reviewed on 12/05/2016

Page 1/10

Ider	ntification
Prod	uct identifier
Trad	e name: Direct-zol RNA PreWash
	<i>le number:</i> R2050-2-40, R2050-2-160 <i>ication of the substance / the mixture</i> Laboratory Reagent
Man Zyme 1706 Irvin USA Phon	<i>ils of the supplier of the safety data sheet</i> <i>ufacturer/Supplier:</i> o Research Corp. 2 Murphy Ave. e, CA 92614 e: 1-949-679-1190 or 1-888-882-9682 ozymoresearch.com
Eme	<i>mation department:</i> Product safety department <i>rgency telephone number:</i> ng normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190
	1/\`1_,•C•,•
	ard(s) identification
Class	ification of the substance or mixture
	GHS02 Flame
Flam	. Liq. 2 H225 Highly flammable liquid and vapor.
	GHS07
Skin	Irrit. 2 H315 Causes skin irritation.
Eye l	Irrit. 2A H319 Causes serious eye irritation.
GHS Haza	<i>l elements</i> <i>label elements</i> The product is classified and labeled according to the Globally Harmonized System (GHS). <i>ard pictograms</i> GHS02, GHS07 <i>al word</i> Danger
Trad Trad <i>Haza</i> High	ard-determining components of labeling: e Secret 001-2100 e Secret 003-2100 ard statements ly flammable liquid and vapor.
Caus	es skin irritation. es serious eye irritation. autionary statements
Keep	away from heat/sparks/open flames/hot surfaces. No smoking.
	explosion-proof electrical/ventilating/lighting/equipment. r protective gloves / eye protection / face protection.
•• Cal	(Contd. on page



Printing date 12/12/2016

Trade name: Direct-zol RNA PreWash

Reviewed on 12/05/2016

Page 2/10

(Contd. of page 1) Ground/bond container and receiving equipment. Keep container tightly closed. Use only non-sparking tools. Take precautionary measures against static discharge. Wash thoroughly after handling. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Specific treatment (see on this label). If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. In case of fire: Use for extinction: CO2, powder or water spray. Take off contaminated clothing and wash it before reuse. Store in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with local/regional/national/international regulations. · Classification system: · NFPA ratings (scale 0 - 4) Health = 1Fire = 3Reactivity = 0· HMIS-ratings (scale 0 - 4) HEALTH 1 Health = 1FIRE 3 Fire = 3Reactivity = 0REACTIVITY 0 · Other hazards · Results of PBT and vPvB assessment · *PBT*: Not applicable. · *vPvB*: Not applicable. **3** Composition/information on ingredients · Chemical characterization: Mixtures · Description: Mixture of the substances listed below with nonhazardous additions.

· Dangerous components:			
Trade Secret 003-2100	≤20%		
Trade Secret 001-2100	≤70%		

4 First-aid measures

· Description of first aid measures

- · After inhalation: Supply fresh air; consult doctor in case of complaints.
- After skin contact: Generally the product does not irritate the skin.

(Contd. on page 3)

US



Printing date 12/12/2016

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

(Contd. of page 2)

Page 3/10

- After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
- After swallowing: Do not induce vomiting; immediately call for medical help.
- Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
- · Suitable extinguishing agents:
- CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- Advice for firefighters
- · Protective equipment: Wear protective clothing.

6 Accidental release measures

- *Personal precautions, protective equipment and emergency procedures* Wear protective equipment. Keep unprotected persons away.
- *Environmental precautions:* Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

- Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- *Reference to other sections*See Section 7 for information on safe handling.
 See Section 8 for information on personal protection equipment.
 See Section 13 for disposal information.

7 Handling and storage

· Handling:

· Precautions for safe handling

No special precautions are necessary if used correctly. Avoid breathing dust, vapor, mist or gas. Avoid prolonged or reapeated exposure. Use only in a chemical fume hood. Open and handle container with care. Keep ignition sources away.

- *Information about protection against explosions and fires:* Keep ignition sources away - Do not smoke. Protect against electrostatic charges.
- · Conditions for safe storage, including any incompatibilities
- · Storage:
- Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: Keep receptacle tightly sealed.

(Contd. on page 4)

⁻ US



Printing date 12/12/2016

Trade name: Direct-zol RNA PreWash

• *Specific end use(s)* No further relevant information available.

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
- · Components with limit values that require monitoring at the workplace:
- The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment:
- · General protective and hygienic measures: Wash hands before breaks and at the end of work.
- · Breathing equipment: Not required.
- · Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation \cdot *Material of gloves*

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

• Eye protection:



Tightly sealed goggles

9 Physical and chemical properties

- · Information on basic physical and chemical properties
- · General Information
- Appearance:

 Form:
 Color:
 Colorless
 Odor:
 Characteristic
 Odor threshold:
 Not determined.

 pH-value:
 Not determined.

(Contd. on page 5)

US

(Contd. of page 3)

Reviewed on 12/05/2016

Page 4/10



Printing date 12/12/2016

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

	(Contd. of page -
• Change in condition Melting point/Melting range: Boiling point/Boiling range:	Undetermined. Undetermined.
· Flash point:	<30 °C (<86 °F)
· Flammability (solid, gaseous):	Not applicable.
· Ignition temperature:	425 °C (797 °F)
· Decomposition temperature:	Not determined.
· Auto igniting:	Product is not selfigniting.
· Danger of explosion:	Product is not explosive. However, formation of explosive air/vapor mixtures are possible.
· Explosion limits: Lower: Upper:	3.5 Vol % 15.0 Vol %
· Vapor pressure at 20 °C (68 °F):	59 hPa (44 mm Hg)
 Density: Relative density Vapor density Evaporation rate 	Not determined. Not determined. Not determined. Not determined.
· Solubility in / Miscibility with	
Water:	Fully miscible.
\cdot Partition coefficient (n-octanol/wate	er): Not determined.
· Viscosity: Dynamic: Kinematic:	Not determined. Not determined.
· Solvent content:	
Organic solvents: VOC content:	69.0 % 69.0 % 690.0 g/l / 5.76 lb/gl
Solids content: • Other information	20.0 % No further relevant information available.

10 Stability and reactivity

· Reactivity No further relevant information available.

· Chemical stability

- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.

(Contd. on page 6)

Page 5/10

US



Printing date 12/12/2016

Trade name: Direct-zol RNA PreWash

· Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- · Information on toxicological effects
- Acute toxicity:
- · LD/LC50 values that are relevant for classification:

Trade Secret 003-2100

Oral LD50 475 mg/kg (rat)

- Primary irritant effect:
- on the skin: No irritant effect.
- on the eye: No irritating effect.
- Sensitization: No sensitizing effects known.
- Additional toxicological information:
- · Carcinogenic categories
- · IARC (International Agency for Research on Cancer)
- Trade Secret 001-2100
- · NTP (National Toxicology Program)

None of the ingredients is listed.

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

· Toxicity

- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · *Bioaccumulative potential* No further relevant information available.
- · Mobility in soil No further relevant information available.
- Additional ecological information:
- · General notes:

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- · Results of PBT and vPvB assessment
- · *PBT*: Not applicable.
- *vPvB*: Not applicable.
- · Other adverse effects No further relevant information available.

Reviewed on 12/05/2016

(Contd. of page 5)

1

(Contd. on page 7)

Page 6/10



Printing date 12/12/2016

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

(Contd. of page 6)

13 Disposal considerations

· Waste treatment methods

· Recommendation:

Dispose of contents in accordance with local/regional/national, and international recommendations.

· Uncleaned packagings:

· Recommendation:

Dispose of container in acoordance with local/regional/national and international recommendations.

· Recommended cleansing agent: Water, if necessary with cleansing agents.

· UN-Number	
· DOT, IMDG, IATA	UN1170
· UN proper shipping name	
$\cdot DOT$	Ethanol mixture
·IMDG	ETHANOL (ETHYL ALCOHOL) mixture
·IATA	ETHANOL mixture
· Transport hazard class(es)	
·DOT	
FLAMMABLE LIQUID	
3	
· Class	3 Flammable liquids
· Label	3
· IMDG, IATA	
V	
· Class	3 Flammable liquids
· Label	3
· Packing group	
· DOT, IMDG, IATA	III
· Environmental hazards:	Not applicable.
· Special precautions for user	Warning: Flammable liquids
· Danger code (Kemler):	30
· EMS Number:	F-E,S-D

Page 7/10



Page 8/10

Safety Data Sheet acc. to OSHA HCS

Printing date 12/12/2016

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

	(Contd. of page 2
· Stowage Category	A
· Transport in bulk according to Annex II of	c · · · · · · · · · · · · · · · · · · ·
MARPOL73/78 and the IBC Code	Not applicable.
· Transport/Additional information:	
·DOT	
· Quantity limitations	On passenger aircraft/rail: 60 L
	On cargo aircraft only: 220 L
· IMDG	
· Limited quantities (LQ)	5L
· Excepted quantities (EQ)	Code: E1
	Maximum net quantity per inner packaging: 30 ml
	Maximum net quantity per outer packaging: 1000 ml
· UN ''Model Regulation'':	UN 1170 ETHANOL MIXTURE, 3, III

15 Regulatory information

 \cdot Safety, health and environmental regulations/legislation specific for the substance or mixture \cdot Sara

· Section 355 (extremely hazardous substances):

None of the ingredients is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

• TSCA (Toxic Substances Control Act):

All ingredients are listed.

· Proposition 65

· Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

 \cdot EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

Trade Secret 001-2100

A3

(Contd. on page 9)

US



Printing date 12/12/2016

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

	(Contd. of page 8		
NIOSH-Ca (National Institute for Occupational Safety and Health)			
None of the ingredients is listed.			
GHS label elements The product is classified Hazard pictograms GHS02, GHS07 Signal word Danger	and labeled according to the Globally Harmonized System (GHS).		
Hazard-determining components of labeling	:		
Trade Secret 001-2100			
Trade Secret 003-2100			
Hazard statements			
Highly flammable liquid and vapor.			
Causes skin irritation.			
Causes serious eye irritation.			
Precautionary statements			
Keep away from heat/sparks/open flames/hot s	surfaces. No smoking.		
Use explosion-proof electrical/ventilating/light			
Wear protective gloves / eye protection / face			
Ground/bond container and receiving equipme	ent.		
Keep container tightly closed.			
Use only non-sparking tools.			
Take precautionary measures against static dis	charge.		
Wash thoroughly after handling.			
	contaminated clothing. Rinse skin with water/shower.		
	eral minutes. Remove contact lenses, if present and easy to do.		
Continue rinsing.			
Specific treatment (see on this label).			
If skin irritation occurs: Get medical advice/at			
If eye irritation persists: Get medical advice/at			
In case of fire: Use for extinction: CO2, powd			
Take off contaminated clothing and wash it be	fore reuse.		
Store in a well-ventilated place. Keep cool.			
	vith local/regional/national/international regulations.		
Chemical safety assessment: A Chemical Safe	ety Assessment has not been carried out.		

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Department issuing SDS: Zymo Research Corp. Safety Department 17062 Murphy Ave. Irvine, CA 92614 USA Phone: 1-949-679-1190 or 1-888-882-9682
Contact: sds@zymoresearch.com
Date of preparation / last revision 12/12/2016 / -

(Contd. on page 10)

Page 9/10

US



Printing date 12/12/2016

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

(Contd. of page 9)

Page 10/10

· Abbreviations and acronyms: ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road) IMDG: International Maritime Code for Dangerous Goods DOT: US Department of Transportation IATA: International Air Transport Association ACGIH: American Conference of Governmental Industrial Hygienists EINECS: European Inventory of Existing Commercial Chemical Substances ELINCS: European List of Notified Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society) NFPA: National Fire Protection Association (USA) HMIS: Hazardous Materials Identification System (USA) VOC: Volatile Organic Compounds (USA, EU) LC50: Lethal concentration, 50 percent LD50: Lethal dose, 50 percent PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative NIOSH: National Institute for Occupational Safety OSHA: Occupational Safety & Health TLV: Threshold Limit Value PEL: Permissible Exposure Limit **REL:** Recommended Exposure Limit Flam. Liq. 2: Flammable liquids - Category 2 Skin Irrit. 2: Skin corrosion/irritation - Category 2 Eye Irrit. 2A: Serious eye damage/eye irritation - Category 2A