

02/16/2017	Kit Components
Product code	Description
R2070 & R2072	Direct-zol RNA MiniPrep Plus, 50 preps and 200 preps
Components:	
R2050-2-40	Direct-zol RNA PreWash
R1003-3-6	RNA Wash Buffer
E1009-A	DNase I
E1010-1-4	DNA Digestion Buffer
W1001-1	DNase/RNase Free Water



Printing date 02/16/2017

Reviewed on 12/05/2016

Page 1/10

Identifi	cation
Product i	dentifier
Trade na	me: Direct-zol RNA PreWash
	umber: R2050-2-40, R2050-2-160 on of the substance / the mixture Laboratory Reagent
Manufact Zymo Res 17062 Mu Irvine, CA USA Phone: 1-	<i>T the supplier of the safety data sheet</i> <i>turer/Supplier:</i> search Corp. urphy Ave. A 92614 949-679-1190 or 1-888-882-9682 oresearch.com
Emergen	<i>ion department:</i> Product safety department <i>cy telephone number:</i> ormal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190
Hazard	(s) identification
Classifica	tion of the substance or mixture
	GHS02 Flame
Flam. Liq	. 2 H225 Highly flammable liquid and vapor.
	GHS07
Skin Irrit.	2 H315 Causes skin irritation.
Eye Irrit.	2A H319 Causes serious eye irritation.
Hazard p	ments el elements The product is classified and labeled according to the Globally Harmonized System (GHS). ictograms GHS02, GHS07 ord Danger
Trade Sec Trade Sec <i>Hazard st</i> Highly fla	etermining components of labeling: cret 001-2100 cret 003-2100 catements ummable liquid and vapor. in irritation.
Causes se <i>Precautio</i>	rious eye irritation. <i>nary statements</i> y from heat/sparks/open flames/hot surfaces. No smoking.
Use explo	sion-proof electrical/ventilating/lighting/equipment.
Wear pro	tective gloves / eye protection / face protection.



Printing date 02/16/2017

Trade name: Direct-zol RNA PreWash

Reviewed on 12/05/2016

(Contd. of page 1) Wear protective gloves / eye protection / face protection. Ground/bond container and receiving equipment. Keep container tightly closed. Use only non-sparking tools. Take precautionary measures against static discharge. Wash thoroughly after handling. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Specific treatment (see on this label). If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. In case of fire: Use for extinction: CO2, powder or water spray. Take off contaminated clothing and wash it before reuse. Store in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with local/regional/national/international regulations. · Classification system: · NFPA ratings (scale 0 - 4) Health = 1Fire = 3Reactivity = 0· HMIS-ratings (scale 0 - 4) HEALTH 1 Health = 1FIRE Fire = 33 Reactivity = 0REACTIVITY 0 · Other hazards · Results of PBT and vPvB assessment · *PBT*: Not applicable. · vPvB: Not applicable. **3** Composition/information on ingredients · Chemical characterization: Mixtures · Description: Mixture of the substances listed below with nonhazardous additions. · Dangerous components:

Dungerous components		
Trade Secret 003-2100	≤20%	
Trade Secret 001-2100	≤70%	

## 4 First-aid measures

· Description of first aid measures

· After inhalation: Supply fresh air; consult doctor in case of complaints.

(Contd. on page 3)

Page 2/10



*Printing date 02/16/2017* 

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

(Contd. of page 2)

Page 3/10

- After skin contact: Generally the product does not irritate the skin.
- · After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
- · After swallowing: Do not induce vomiting; immediately call for medical help.
- · Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

#### 5 Fire-fighting measures

- · Extinguishing media
- Suitable extinguishing agents:
- CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- Advice for firefighters
- · Protective equipment: Wear protective clothing.

## 6 Accidental release measures

- *Personal precautions, protective equipment and emergency procedures* Wear protective equipment. Keep unprotected persons away.
- Environmental precautions: Dilute with plenty of water.
- Do not allow to enter sewers/ surface or ground water.
- Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- · Reference to other sections
- See Section 7 for information on safe handling.
- See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.

## 7 Handling and storage

- · Handling:
- · Precautions for safe handling

No special precautions are necessary if used correctly. Avoid breathing dust, vapor, mist or gas. Avoid prolonged or reapeated exposure. Use only in a chemical fume hood. Open and handle container with care. Keep ignition sources away.

- *Information about protection against explosions and fires:* Keep ignition sources away - Do not smoke.
- Protect against electrostatic charges.
- Conditions for safe storage, including any incompatibilities • Storage:
- Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: Keep receptacle tightly sealed.

(Contd. on page 4)



*Printing date 02/16/2017* 

#### Trade name: Direct-zol RNA PreWash

• *Specific end use*(*s*) No further relevant information available.

## 8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
- · Components with limit values that require monitoring at the workplace:
- The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment:
- · General protective and hygienic measures: Wash hands before breaks and at the end of work.
- · Breathing equipment: Not required.
- · Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation  $\cdot$  *Material of gloves* 

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

#### · Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

• Eye protection:



Tightly sealed goggles

## 9 Physical and chemical properties

- · Information on basic physical and chemical properties
- · General Information

• Appearance:	
Form:	Fluid
Color:	Colorless
· Odor:	Characteristic
• Odor threshold:	Not determined.
· pH-value:	Not determined.

(Contd. on page 5)

Page 4/10

(Contd. of page 3)

Reviewed on 12/05/2016

US



Printing date 02/16/2017

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

	(Contd. of page
Change in condition	
Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	Undetermined.
Flash point:	<30 °C (<86 °F)
Flammability (solid, gaseous):	Not applicable.
Ignition temperature:	425 °C (797 °F)
Decomposition temperature:	Not determined.
Auto igniting:	Product is not selfigniting.
Danger of explosion:	Product is not explosive. However, formation of explosive air/vapor mixtures are possible.
Explosion limits:	
Lower:	3.5 Vol %
Upper:	15.0 Vol %
Vapor pressure at 20 °C (68 °F):	59 hPa (44 mm Hg)
Density:	Not determined.
Relative density	Not determined.
Vapor density	Not determined.
Evaporation rate	Not determined.
Solubility in / Miscibility with	
Water:	Fully miscible.
Partition coefficient (n-octanol/wate	er): Not determined.
Viscosity:	
Dynamic:	Not determined.
Kinematic:	Not determined.
Solvent content:	
Organic solvents:	69.0 %
VOC content:	69.0 %
	690.0 g/l / 5.76 lb/gl
Solids content:	20.0 %
Other information	No further relevant information available.

# 10 Stability and reactivity

· Reactivity No further relevant information available.

· Chemical stability

- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.

(Contd. on page 6)

Page 5/10

US



Printing date 02/16/2017

Trade name: Direct-zol RNA PreWash

· Hazardous decomposition products: No dangerous decomposition products known.

## 11 Toxicological information

- · Information on toxicological effects
- Acute toxicity:
- · LD/LC50 values that are relevant for classification:

Trade Secret 003-2100

Oral LD50 475 mg/kg (rat)

- · Primary irritant effect:
- · on the skin: No irritant effect.
- · on the eye: No irritating effect.
- · Sensitization: No sensitizing effects known.
- · Additional toxicological information:

· Carcinogenic categories

· IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

# 12 Ecological information

· Toxicity

- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes:

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- · Results of PBT and vPvB assessment
- · PBT: Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

(Contd. on page 7)

Reviewed on 12/05/2016

(Contd. of page 5)

Page 6/10



Printing date 02/16/2017

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

(Contd. of page 6)

# 13 Disposal considerations

· Waste treatment methods

· Recommendation:

Dispose of contents in accordance with local/regional/national, and international recommendations.

· Uncleaned packagings:

· Recommendation:

Dispose of container in acoordance with local/regional/national and international recommendations.

· Recommended cleansing agent: Water, if necessary with cleansing agents.

· UN-Number	
· DOT, IMDG, IATA	UN1170
· UN proper shipping name	
$\cdot DOT$	Ethanol mixture
·IMDG	ETHANOL (ETHYL ALCOHOL) mixture
·IATA	ETHANOL mixture
· Transport hazard class(es)	
·DOT	
· Class	3 Flammable liquids
· Label	3
· IMDG, IATA	
3	
· Class	3 Flammable liquids
· Label	3
· Packing group	
· DOT, IMDG, IATA	III
• Environmental hazards:	Not applicable.
$\cdot$ Special precautions for user	Warning: Flammable liquids
· Danger code (Kemler): · EMS Number:	30 F-E,S-D



Page 8/10

# Safety Data Sheet acc. to OSHA HCS

Printing date 02/16/2017

Reviewed on 12/05/2016

#### Trade name: Direct-zol RNA PreWash

	(Contd. of pag
· Stowage Category	А
• Transport in bulk according to Annex	II of
MARPOL73/78 and the IBC Code	Not applicable.
· Transport/Additional information:	
DOT	
Quantity limitations	On passenger aircraft/rail: 60 L
	On cargo aircraft only: 220 L
· IMDG	
Limited quantities (LQ)	5L
Excepted quantities (EQ)	Code: E1
	Maximum net quantity per inner packaging: 30 ml
	Maximum net quantity per outer packaging: 1000 ml
· UN ''Model Regulation'':	UN 1170 ETHANOL MIXTURE, 3, III

# 15 Regulatory information

 $\cdot$  Safety, health and environmental regulations/legislation specific for the substance or mixture  $\cdot$  Sara

· Section 355 (extremely hazardous substances):

None of the ingredients is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

• TSCA (Toxic Substances Control Act):

All ingredients are listed.

· Proposition 65

· Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

 $\cdot$  EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

Trade Secret 001-2100

A3

(Contd. on page 9)

US



Printing date 02/16/2017

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

	(Contd. of page 8
NIOSH-Ca (National Institute for Occupational Safety and Health)	
None of the ingredients is listed.	
GHS label elements The product is classified and labeled according to Hazard pictograms GHS02, GHS07 Signal word Danger	the Globally Harmonized System (GHS).
Hazard-determining components of labeling:	
Trade Secret 001-2100	
Trade Secret 003-2100	
Hazard statements	
Highly flammable liquid and vapor.	
Causes skin irritation.	
Causes serious eye irritation.	
Precautionary statements	
Keep away from heat/sparks/open flames/hot surfaces. No smoking.	
Use explosion-proof electrical/ventilating/lighting/equipment.	
Wear protective gloves / eye protection / face protection.	
Wear protective gloves / eye protection / face protection.	
Ground/bond container and receiving equipment.	
Keep container tightly closed.	
Use only non-sparking tools.	
Take precautionary measures against static discharge. Wash thoroughly after handling.	
	ince alin with water/abower
If on skin (or hair): Take off immediately all contaminated clothing. Ri If in eyes: Rinse cautiously with water for several minutes. Remove con	
Continue rinsing.	mact lenses, if present and easy to do.
Specific treatment (see on this label).	
If skin irritation occurs: Get medical advice/attention.	
If eye irritation persists: Get medical advice/attention.	
In case of fire: Use for extinction: CO2, powder or water spray.	
Take off contaminated clothing and wash it before reuse.	
Store in a well-ventilated place. Keep cool.	
Dispose of contents/container in accordance with local/regional/nation	al/international regulations
<b>Chemical safety assessment:</b> A Chemical Safety Assessment has not b	

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

• Department issuing SDS:

Zymo Research Corp. Safety Department 17062 Murphy Ave. Irvine, CA 92614 USA Phone: 1-949-679-1190 or 1-888-882-9682 · *Contact:* sds@zymoresearch.com

(Contd. on page 10)

Page 9/10

US



Printing date 02/16/2017

Reviewed on 12/05/2016

Trade name: Direct-zol RNA PreWash

(Contd. of page 9) · Date of preparation / last revision 02/16/2017 / -• Abbreviations and acronyms: ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road) IMDG: International Maritime Code for Dangerous Goods DOT: US Department of Transportation IATA: International Air Transport Association ACGIH: American Conference of Governmental Industrial Hygienists EINECS: European Inventory of Existing Commercial Chemical Substances ELINCS: European List of Notified Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society) NFPA: National Fire Protection Association (USA) HMIS: Hazardous Materials Identification System (USA) VOC: Volatile Organic Compounds (USA, EU) LC50: Lethal concentration, 50 percent LD50: Lethal dose, 50 percent PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative NIOSH: National Institute for Occupational Safety OSHA: Occupational Safety & Health TLV: Threshold Limit Value PEL: Permissible Exposure Limit REL: Recommended Exposure Limit Flam. Liq. 2: Flammable liquids - Category 2 Skin Irrit. 2: Skin corrosion/irritation - Category 2 Eye Irrit. 2A: Serious eye damage/eye irritation - Category 2A

Page 10/10



Printing date 02/16/2017

Reviewed on 01/22/2016

inting date 02/16/2017	Reviewed on 01/22/2016
I Identification	
Product identifier	
Trade name: RNA Wash Buffer	
Article number: R1003-3-6, R1003-3-12, R1003-3-24, R100 Application of the substance / the mixture Laboratory Reag	
Details of the supplier of the safety data sheet Manufacturer/Supplier: Zymo Research Corp. 17062 Murphy Ave. Irvine, CA 92614 USA Phone: 1-949-679-1190 or 1-888-882-9682 sds@zymoresearch.com	
<ul> <li><i>Information department:</i> Product safety department</li> <li><i>Emergency telephone number:</i> During normal business hours (8 am to 5 pm Pacific Standard)</li> </ul>	d Time): +1 (949) 679 1190
· Classification of the substance or mixture	onized System (GHS).
<ul> <li><i>Pazard(s) identification</i></li> <li><i>Classification of the substance or mixture</i> The product is not classified according to the Globally Harm</li> <li><i>Label elements</i></li> <li><i>GHS label elements</i> Void</li> <li><i>Hazard pictograms</i> Void</li> </ul>	onized System (GHS).
<ul> <li>Classification of the substance or mixture The product is not classified according to the Globally Harm</li> <li>Label elements</li> <li>GHS label elements Void</li> <li>Hazard pictograms Void</li> <li>Signal word Void</li> </ul>	ionized System (GHS).
<ul> <li>Classification of the substance or mixture The product is not classified according to the Globally Harm</li> <li>Label elements</li> <li>GHS label elements Void</li> <li>Hazard pictograms Void</li> <li>Signal word Void</li> <li>Hazard statements Void</li> <li>Classification system:</li> </ul>	ionized System (GHS).
Classification of the substance or mixture The product is not classified according to the Globally Harm Label elements GHS label elements Void Hazard pictograms Void Signal word Void Hazard statements Void Classification system:	ionized System (GHS).
<ul> <li>Classification of the substance or mixture The product is not classified according to the Globally Harm</li> <li>Label elements</li> <li>GHS label elements Void</li> <li>Hazard pictograms Void</li> <li>Signal word Void</li> <li>Hazard statements Void</li> <li>Classification system:</li> <li>NFPA ratings (scale 0 - 4)</li> <li>Health = 0 Fire = 0</li> </ul>	ionized System (GHS).
Classification of the substance or mixture The product is not classified according to the Globally Harm Label elements GHS label elements Void Hazard pictograms Void Signal word Void Hazard statements Void Classification system: NFPA ratings (scale $0 - 4$ ) Health = 0 Fire = 0 Reactivity = 0	onized System (GHS).
Classification of the substance or mixture The product is not classified according to the Globally Harm Label elements GHS label elements Void Hazard pictograms Void Signal word Void Hazard statements Void Classification system: NFPA ratings (scale $0 - 4$ ) Health = 0 Fire = 0 Reactivity = 0 HMIS-ratings (scale $0 - 4$ ) Health = 0 Fire = 0 Reactivity = 0 Health = 0 Fire = 0 Fire = 0 Health = 0 Fire = 0	onized System (GHS).

(Contd. on page 2)

Page 1/8



Printing date 02/16/2017

Trade name: RNA Wash Buffer

Reviewed on 01/22/2016

(Contd. of page 1)

Page 2/8

## 3 Composition/information on ingredients

- · Chemical characterization: Mixtures
- · Description: Mixture of the substances listed below with nonhazardous additions.
- · Dangerous components: Void

#### 4 First-aid measures

- · Description of first aid measures
- · General information: No special measures required.
- · After inhalation: Supply fresh air; consult doctor in case of complaints.
- · After skin contact: Generally the product does not irritate the skin.
- · After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
- · After swallowing: Do not induce vomiting; immediately call for medical help.
- Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

## 5 Fire-fighting measures

- · Extinguishing media
- · Suitable extinguishing agents:
- CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- · Advice for firefighters
- · Protective equipment: Wear protective clothing.

## 6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Wear protective clothing.
- · Environmental precautions:
- Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

- · Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- · Reference to other sections
- See Section 7 for information on safe handling.
- See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.
- · Protective Action Criteria for Chemicals

#### • PAC-1:

CAS: 1185-53-1	2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride	12 mg/m3
CAS: 6381-92-6	Edetate Disodium, Dihydrate	30 mg/m3
	(Co	ntd. on page 3)



Printing date 02/16/2017

Reviewed on 01/22/2016

#### Trade name: RNA Wash Buffer

		(Contd. of page 2)
· PAC-2:		
CAS: 1185-53-1	2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride	130 mg/m3
CAS: 6381-92-6	Edetate Disodium, Dihydrate	330 mg/m3
· PAC-3:		
CAS: 1185-53-1	2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride	790 mg/m3
CAS: 6381-92-6	Edetate Disodium, Dihydrate	2,000 mg/m3

## 7 Handling and storage

· Handling:

- · Precautions for safe handling No special measures required.
- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities

· Storage:

- Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: None.
- · Specific end use(s) Laboratory reagent

## 8 Exposure controls/personal protection

- · Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
- · Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

• Additional information: The lists that were valid during the creation were used as basis.

- · Exposure controls
- · Personal protective equipment:
- General protective and hygienic measures:
- The usual precautionary measures for handling chemicals should be followed.
- · Breathing equipment: Not required.
- Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation  $\cdot$  *Material of gloves* 

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

(Contd. on page 4)

Page 3/8

US



Printing date 02/16/2017

Reviewed on 01/22/2016

(Contd. of page 3)

## Trade name: RNA Wash Buffer

#### · Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

• *Eye protection:* Goggles recommended during refilling.

Information on basic physical and	chemical properties	
General Information Appearance:		
Appearance: Form:	Liquid	
Color:	Clear	
Odor:	Odorless	
Odor threshold:	Not determined.	
pH-value:	Not determined.	
Change in condition		
Melting point/Melting range:	Undetermined.	
Boiling point/Boiling range:	Undetermined.	
Flash point:	Not applicable.	
Flammability (solid, gaseous):	Not applicable.	
Ignition temperature:		
Decomposition temperature:	Not determined.	
Auto igniting:	Product is not selfigniting.	
Danger of explosion:	Product does not present an explosion hazard.	
Explosion limits:		
Lower:	Not Applicable	
Upper:	Not Applicable	
Vapor pressure:	Not determined.	
Density:	Not determined.	
Relative density	Not determined.	
Vapor density	Not determined.	
Evaporation rate	Not determined.	
Solubility in / Miscibility with		
Water:	Fully miscible.	
Partition coefficient (n-octanol/wat	er): Not determined.	
Viscosity:		
Dynamic:	Not determined.	
Kinematic:	Not determined.	



Printing date 02/16/2017

Reviewed on 01/22/2016

n

## 10 Stability and reactivity

- · Reactivity No further relevant information available.
- · Chemical stability
- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

## 11 Toxicological information

- · Information on toxicological effects
- · Acute toxicity:
- · Primary irritant effect:
- on the skin: No irritant effect.
- on the eye: No irritating effect.
- · Sensitization: No sensitizing effects known.
- Additional toxicological information:

The product is not subject to classification according to internally approved calculation methods for preparations: When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

#### · Carcinogenic categories

· IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

- · NTP (National Toxicology Program)
- None of the ingredients is listed.
- · OSHA-Ca (Occupational Safety & Health Administration)
- None of the ingredients is listed.

## 12 Ecological information

- · Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.

(Contd. on page 6)

Page 5/8



*Printing date 02/16/2017* 

Reviewed on 01/22/2016

(Contd. of page 5)

#### Trade name: RNA Wash Buffer

· Behavior in environmental systems:

- · *Bioaccumulative potential* No further relevant information available.
- · Mobility in soil No further relevant information available.
- Additional ecological information:
- · General notes:

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- · Results of PBT and vPvB assessment
- *PBT*: Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

## 13 Disposal considerations

- · Waste treatment methods
- · *Recommendation:* Smaller quantities can be disposed of with household waste.
- · Uncleaned packagings:
- · Recommendation:
- Dispose of container in acoordance with local/regional/national and international recommendations.
- · Recommended cleansing agent: Water, if necessary with cleansing agents.

UN-Number DOT, ADN, IMDG, IATA	not regulated
UN proper shipping name DOT, ADN, IMDG, IATA	not regulated
Transport hazard class(es)	
DOT, ADN, IMDG, IATA Class	not regulated
Packing group DOT, IMDG, IATA	not regulated
Environmental hazards:	Not applicable.
Special precautions for user	Not applicable.
Transport in bulk according to Annex MARPOL73/78 and the IBC Code	II of Not applicable.
UN ''Model Regulation'':	not regulated

(Contd. on page 7)

Page 6/8



Printing date 02/16/2017

Trade name: RNA Wash Buffer

Reviewed on 01/22/2016

(Contd. of page 6)

Page 7/8

## **15 Regulatory information**

 $\cdot$  Safety, health and environmental regulations/legislation specific for the substance or mixture  $\cdot$  Sara

· Section 355 (extremely hazardous substances):

None of the ingredients is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

· TSCA (Toxic Substances Control Act):

CAS: 1185-53-1 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride

· Proposition 65

· Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

· GHS label elements Void

· Hazard pictograms Void

· Signal word Void

· Hazard statements Void

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

# 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS:

Zymo Research Corp. Safety Department 17062 Murphy Ave. Irvine, CA 92614 USA



Printing date 02/16/2017

Reviewed on 01/22/2016


#### Trade name: RNA Wash Buffer

(Contd. of page 7) Phone: 1-949-679-1190 or 1-888-882-9682 · Contact: sds@zymoresearch.com · Date of preparation / last revision 02/16/2017 / -· Abbreviations and acronyms: ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road) IMDG: International Maritime Code for Dangerous Goods DOT: US Department of Transportation IATA: International Air Transport Association ACGIH: American Conference of Governmental Industrial Hygienists EINECS: European Inventory of Existing Commercial Chemical Substances ELINCS: European List of Notified Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society) NFPA: National Fire Protection Association (USA) HMIS: Hazardous Materials Identification System (USA) VOC: Volatile Organic Compounds (USA, EU) PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative NIOSH: National Institute for Occupational Safety OSHA: Occupational Safety & Health TLV: Threshold Limit Value PEL: Permissible Exposure Limit REL: Recommended Exposure Limit



*Printing date 02/16/2017* 

Reviewed on 02/02/2016



Page 1/8



*Printing date 02/16/2017* 

Reviewed on 02/02/2016

Trade name: DNase I

(Contd. of page 1)

Page 2/8

#### 3 Composition/information on ingredients

- · Chemical characterization: Mixtures
- · Description: Mixture of the substances listed below with nonhazardous additions.
- · Dangerous components: Void

#### 4 First-aid measures

- · Description of first aid measures
- · General information: No special measures required.
- · After inhalation: Supply fresh air; consult doctor in case of complaints.
- · After skin contact: Generally the product does not irritate the skin.
- · After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
- · After swallowing: Do not induce vomiting; immediately call for medical help.
- Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

## 5 Fire-fighting measures

- · Extinguishing media
- · Suitable extinguishing agents:
- CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- · Advice for firefighters
- · Protective equipment: Wear protective clothing.

## 6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Wear protective clothing.
- · Environmental precautions: Dilute with plenty of water.
- · Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- · Reference to other sections
- See Section 7 for information on safe handling.
- See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.
- · Protective Action Criteria for Chemicals

· PAC-1:

None of the ingredients is listed.

· PAC-2:

None of the ingredients is listed.

(Contd. on page 3)



*Printing date 02/16/2017* 

Reviewed on 02/02/2016

(Contd. of page 2)

#### Trade name: DNase I

#### • PAC-3:

None of the ingredients is listed.

## 7 Handling and storage

#### · Handling:

- · *Precautions for safe handling* No special measures required.
- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- · Storage:
- Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Not required.
- Further information about storage conditions: None.
- Specific end use(s) Laboratory reagent

## 8 Exposure controls/personal protection

• Additional information about design of technical systems: No further data; see item 7.

- · Control parameters
- · Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- · Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment:
- · General protective and hygienic measures:

The usual precautionary measures for handling chemicals should be followed.

- · Breathing equipment: Not required.
- Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation  $\cdot$  *Material of gloves* 

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

• Eye protection: Goggles recommended during refilling.

(Contd. on page 4)

Page 3/8



Printing date 02/16/2017

Trade name: DNase I

· Appearance: Form:

Color:

· pH-value:

· Flash point:

• Auto igniting:

· Odor:

#### 9 Physical and chemical properties · Information on basic physical and chemical properties · General Information Solid White Odorless · Odor threshold: Not determined. Not determined. · Change in condition Melting point/Melting range: Undetermined. Boiling point/Boiling range: Undetermined. Not applicable. · Flammability (solid, gaseous): Not applicable. · Ignition temperature: **Decomposition temperature:** Not determined. Product is not selfigniting. · Danger of explosion: Product does not present an explosion hazard. · Explosion limits: Not Applicable Not Applicable

Lower: Upper: · Vapor pressure: Not determined. · Density: Not determined. · Relative density Not determined. · Vapor density Not determined. · Evaporation rate Not determined. · Solubility in / Miscibility with Water: Fully miscible. · Partition coefficient (n-octanol/water): Not determined. · Viscosity: Dynamic: Not determined. Kinematic: Not determined. · Solvent content: Organic solvents: 0.0 % **VOC content:** 0.0 g/l / 0.00 lb/gl · Other information No further relevant information available.



Reviewed on 02/02/2016

(Contd. of page 3)

(Contd. on page 5)

US



*Printing date 02/16/2017* 

Reviewed on 02/02/2016

#### Trade name: DNase I

#### (Contd. of page 4)

Page 5/8

#### 10 Stability and reactivity

- · *Reactivity* No further relevant information available.
- · Chemical stability
- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

## 11 Toxicological information

- · Information on toxicological effects
- Acute toxicity:
- · Primary irritant effect:
- on the skin: No irritant effect.
- on the eye: No irritating effect.
- Sensitization: No sensitizing effects known.
- · Additional toxicological information:

The product is not subject to classification according to internally approved calculation methods for preparations: When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

#### · Carcinogenic categories

· IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

#### · OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

## 12 Ecological information

- · Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes: Not known to be hazardous to water.
- · Results of PBT and vPvB assessment
- · *PBT*: Not applicable.
- · *vPvB*: Not applicable.

(Contd. on page 6)



*Printing date 02/16/2017* 

Reviewed on 02/02/2016

(Contd. of page 5)

Page 6/8

#### Trade name: DNase I

· Other adverse effects No further relevant information available.

#### 13 Disposal considerations

- · Waste treatment methods
- · *Recommendation:* Smaller quantities can be disposed of with household waste.
- · Uncleaned packagings:
- · Recommendation:
- Dispose of container in acoordance with local/regional/national and international recommendations.
- · Recommended cleansing agent: Water, if necessary with cleansing agents.

UN-Number DOT, ADN, IMDG, IATA	not regulated
UN proper shipping name DOT, ADN, IMDG, IATA	not regulated
Transport hazard class(es)	
DOT, ADN, IMDG, IATA Class	not regulated
Packing group DOT, IMDG, IATA	not regulated
Environmental hazards:	Not applicable.
Special precautions for user	Not applicable.
Transport in bulk according to Annex I MARPOL73/78 and the IBC Code	II of Not applicable.
UN "Model Regulation":	not regulated

# 15 Regulatory information

 $\cdot$  Safety, health and environmental regulations/legislation specific for the substance or mixture  $\cdot$  Sara

· Section 355 (extremely hazardous substances):

None of the ingredients is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

• *TSCA (Toxic Substances Control Act):* None of the ingredients is listed.

(Contd. on page 7)

US



Printing date 02/16/2017

Reviewed on 02/02/2016

Trade name: DNase I

· Proposition 65

(Contd. of page 6)

Page 7/8

· Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

· GHS label elements Void

· Hazard pictograms Void

· Signal word Void

· Hazard statements Void

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS:

Zymo Research Corp. Safety Department 17062 Murphy Ave. Irvine, CA 92614 USA Phone: 1-949-679-1190 or 1-888-882-9682 · Contact: sds@zymoresearch.com · Date of preparation / last revision 02/16/2017 / -· Abbreviations and acronyms: ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road) IMDG: International Maritime Code for Dangerous Goods DOT: US Department of Transportation IATA: International Air Transport Association ACGIH: American Conference of Governmental Industrial Hygienists EINECS: European Inventory of Existing Commercial Chemical Substances ELINCS: European List of Notified Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society)

(Contd. on page 8)

US



*Printing date 02/16/2017* 

#### Trade name: DNase I

Reviewed on 02/02/2016

NFPA: National Fire Protection Association (USA) HMIS: Hazardous Materials Identification System (USA) VOC: Volatile Organic Compounds (USA, EU) PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative NIOSH: National Institute for Occupational Safety OSHA: Occupational Safety & Health TLV: Threshold Limit Value PEL: Permissible Exposure Limit REL: Recommended Exposure Limit (Contd. of page 7)

us –


Page 8/8



Printing date 02/16/2017

Reviewed on 02/03/2016

Page 1/8





Printing date 02/16/2017

Trade name: DNA Digestion Buffer

Reviewed on 02/03/2016

(Contd. of page 1)

## 3 Composition/information on ingredients

- · Chemical characterization: Mixtures
- · Description: Mixture of the substances listed below with nonhazardous additions.
- · Dangerous components: Void

#### 4 First-aid measures

- · Description of first aid measures
- · General information: No special measures required.
- · After inhalation: Supply fresh air; consult doctor in case of complaints.
- · After skin contact: Generally the product does not irritate the skin.
- · After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
- After swallowing: Do not induce vomiting; immediately call for medical help.
- Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

## 5 Fire-fighting measures

- · Extinguishing media
- Suitable extinguishing agents:
- CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- · Advice for firefighters
- · Protective equipment: Wear protective clothing.

## 6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Wear protective clothing.
- · Environmental precautions: Dilute with plenty of water.
- · Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- · Reference to other sections
- See Section 7 for information on safe handling.
- See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.
- · Protective Action Criteria for Chemicals

· PAC-1:

None of the ingredients is listed.

· PAC-2:

None of the ingredients is listed.

(Contd. on page 3)

#### Page 2/8



Printing date 02/16/2017

Reviewed on 02/03/2016

(Contd. of page 2)

#### Trade name: DNA Digestion Buffer

#### · PAC-3:

None of the ingredients is listed.

## 7 Handling and storage

#### · Handling:

- · *Precautions for safe handling* No special measures required.
- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- · Storage:
- Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Not required.
- Further information about storage conditions: None.
- Specific end use(s) Laboratory reagent

## 8 Exposure controls/personal protection

• Additional information about design of technical systems: No further data; see item 7.

- · Control parameters
- · Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- · Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment:
- · General protective and hygienic measures:

The usual precautionary measures for handling chemicals should be followed.

- · Breathing equipment: Not required.
- Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation  $\cdot$  *Material of gloves* 

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

• Eye protection: Goggles recommended during refilling.

(Contd. on page 4)

Page 3/8



Printing date 02/16/2017

Trade name: DNA Digestion Buffer

Reviewed on 02/03/2016

(Contd. of page 3)

Information on basic physical and	chemical properties
General Information	enemieu properaes
Appearance:	
Form:	Liquid
Color:	Colorless
Odor:	Odorless
Odor threshold:	Not determined.
pH-value:	Not determined.
Change in condition	
Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	Undetermined.
Flash point:	Not applicable.
Flammability (solid, gaseous):	Not applicable.
Ignition temperature:	
Decomposition temperature:	Not determined.
Auto igniting:	Product is not selfigniting.
Danger of explosion:	Product does not present an explosion hazard.
Explosion limits:	
Lower:	Not Applicable
Upper:	Not Applicable
Vapor pressure:	Not determined.
Density:	Not determined.
Relative density	Not determined.
Vapor density	Not determined.
Evaporation rate	Not determined.
Solubility in / Miscibility with	
Water:	Fully miscible.
Partition coefficient (n-octanol/wat	ter): Not determined.
Viscosity:	
Dynamic:	Not determined.
Kinematic:	Not determined.
Solvent content:	
Organic solvents:	0.0 %
VOC content:	0.0 g/l / 0.00 lb/gl
Other information	No further relevant information available.

(Contd. on page 5)

US

*Page 4/8* 



Printing date 02/16/2017

Reviewed on 02/03/2016

Trade name: DNA Digestion Buffer

(Contd. of page 4)

Page 5/8

## 10 Stability and reactivity

- · *Reactivity* No further relevant information available.
- · Chemical stability
- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

## 11 Toxicological information

- · Information on toxicological effects
- Acute toxicity:
- · Primary irritant effect:
- on the skin: No irritant effect.
- on the eye: No irritating effect.
- Sensitization: No sensitizing effects known.
- · Additional toxicological information:

The product is not subject to classification according to internally approved calculation methods for preparations: When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

#### · Carcinogenic categories

· IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

#### · OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

# 12 Ecological information

- · Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes: Generally not hazardous for water
- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- *vPvB*: Not applicable.

(Contd. on page 6)

US



*Printing date 02/16/2017* 

Trade name: DNA Digestion Buffer

· Other adverse effects No further relevant information available.

## 13 Disposal considerations

· Waste treatment methods

· *Recommendation:* Smaller quantities can be disposed of with household waste.

· Uncleaned packagings:

· Recommendation:

Dispose of container in acoordance with local/regional/national and international recommendations.

· Recommended cleansing agent: Water, if necessary with cleansing agents.

UN-Number	
DOT, ADN, IMDG, IATA	not regulated
UN proper shipping name DOT, ADN, IMDG, IATA	not regulated
Transport hazard class(es)	
DOT, ADN, IMDG, IATA	
Class	not regulated
Packing group	
DOT, IMDG, IATA	not regulated
Environmental hazards:	Not applicable.
Special precautions for user	Not applicable.
Transport in bulk according to Annex II	l of
MARPOL73/78 and the IBC Code	Not applicable.

# 15 Regulatory information

 $\cdot$  Safety, health and environmental regulations/legislation specific for the substance or mixture  $\cdot$  Sara

· Section 355 (extremely hazardous substances):

None of the ingredients is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

• *TSCA (Toxic Substances Control Act):* None of the ingredients is listed.

(Contd. on page 7)

Page 6/8

*Reviewed on 02/03/2016* 

(Contd. of page 5)

US



Printing date 02/16/2017

Reviewed on 02/03/2016

Trade name: DNA Digestion Buffer

(Contd. of page 6)

Page 7/8

Proposition 65
 Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

· GHS label elements Void

· Hazard pictograms Void

· Signal word Void

· Hazard statements Void

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS:

Zymo Research Corp. Safety Department 17062 Murphy Ave. Irvine, CA 92614 USA Phone: 1-949-679-1190 or 1-888-882-9682 · Contact: sds@zymoresearch.com · Date of preparation / last revision 02/16/2017 / -· Abbreviations and acronyms: ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road) IMDG: International Maritime Code for Dangerous Goods DOT: US Department of Transportation IATA: International Air Transport Association ACGIH: American Conference of Governmental Industrial Hygienists EINECS: European Inventory of Existing Commercial Chemical Substances ELINCS: European List of Notified Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society)

(Contd. on page 8)

US



Printing date 02/16/2017

Trade name: DNA Digestion Buffer

Reviewed on 02/03/2016

NFPA: National Fire Protection Association (USA) HMIS: Hazardous Materials Identification System (USA) VOC: Volatile Organic Compounds (USA, EU) PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative NIOSH: National Institute for Occupational Safety OSHA: Occupational Safety & Health TLV: Threshold Limit Value PEL: Permissible Exposure Limit REL: Recommended Exposure Limit (Contd. of page 7)

US –

Page 8/8



Printing date 02/16/2017

Reviewed on 12/03/2015

Identification	
Product identifier	
· Trade name: DNase/RNase Free Water	
<ul> <li>Article number: W1001-1, W1001-4, W1001-6, W1001-10, W1001-30</li> <li>CAS Number: 7732-18-5</li> <li>EC number: 231-791-2</li> <li>Application of the substance / the mixture Laboratory Reagent</li> </ul>	
<ul> <li>Details of the supplier of the safety data sheet</li> <li>Manufacturer/Supplier: Zymo Research Corp.</li> <li>17062 Murphy Ave.</li> <li>Irvine, CA 92614</li> <li>USA</li> <li>Phone: 1-949-679-1190 or 1-888-882-9682</li> <li>sds@zymoresearch.com</li> </ul>	
<ul> <li>Information department: Product safety department</li> <li>Emergency telephone number: During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190</li> </ul>	
<ul> <li><i>Plazard(s) identification</i></li> <li><i>Classification of the substance or mixture</i> The substance is not classified according to the Globally Harmonized System (GHS).</li> </ul>	
<ul> <li>Classification of the substance or mixture The substance is not classified according to the Globally Harmonized System (GHS).</li> <li>Label elements</li> <li>GHS label elements Void</li> <li>Hazard pictograms Void</li> <li>Signal word Void</li> <li>Hazard statements Void</li> <li>Classification system:</li> </ul>	
<ul> <li>Hazard(s) identification</li> <li>Classification of the substance or mixture The substance is not classified according to the Globally Harmonized System (GHS).</li> <li>Label elements</li> <li>GHS label elements Void</li> <li>Hazard pictograms Void</li> <li>Signal word Void</li> <li>Hazard statements Void</li> <li>Classification system:</li> </ul>	
<ul> <li><i>Hazard(s) identification</i></li> <li><i>Classification of the substance or mixture</i> The substance is not classified according to the Globally Harmonized System (GHS).</li> <li><i>Label elements</i></li> <li><i>GHS label elements</i> Void</li> <li><i>Hazard pictograms</i> Void</li> <li><i>Signal word</i> Void</li> <li><i>Hazard statements</i> Void</li> <li><i>Classification system:</i></li> <li><i>NFPA ratings (scale 0 - 4)</i></li> <li>Health = 0 Fire = 0</li> </ul>	
$\begin{array}{l} \textbf{Hazard(s) identification} \\ \hline \textbf{Classification of the substance or mixture} \\ \hline \textbf{The substance is not classified according to the Globally Harmonized System (GHS).} \\ \hline \textbf{Label elements} \\ \hline \textbf{GHS label elements Void} \\ \hline \textbf{Hazard pictograms Void} \\ \hline \textbf{Signal word Void} \\ \hline \textbf{Hazard statements Void} \\ \hline \textbf{Classification system:} \\ \hline \textbf{NFPA ratings (scale 0 - 4)} \\ \hline \textbf{Health} = 0 \\ \hline \textbf{Fire} = 0 \\ \hline \textbf{Reactivity} = 0 \\ \end{array}$	
<b>Classification of the substance or mixture</b> The substance is not classified according to the Globally Harmonized System (GHS).         • Label elements         • GHS label elements         • GHS label elements Void         • Hazard pictograms Void         • Signal word Void         • Hazard statements Void         • Classification system:         • NFPA ratings (scale 0 - 4)         • Mis-ratings (scale 0 - 4)         • Health = 0         Fire = 0         Reactivity = 0         • HMIS-ratings (scale 0 - 4)         • Health = 0         Fire = 0         Reactivity = 0	

Page 1/8



Printing date 02/16/2017

Trade name: DNase/RNase Free Water

· vPvB: Not applicable.

## 3 Composition/information on ingredients

- · Chemical characterization: Substances
- · CAS No. Description
- 7732-18-5 water, distilled, conductivity or of similar purity
- · Identification number(s)
- *EC number:* 231-791-2

## 4 First-aid measures

- · Description of first aid measures
- · General information: No special measures required.
- · After inhalation: Supply fresh air; consult doctor in case of complaints.
- After skin contact: Generally the product does not irritate the skin.
- · After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
- After swallowing: Do not induce vomiting; immediately call for medical help.
- Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

## 5 Fire-fighting measures

- · Extinguishing media
- · Suitable extinguishing agents:
- CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam. • *Advice for firefighters*
- · Protective equipment: Wear protective clothing.

## 6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Wear protective clothing.
- · Environmental precautions: Dilute with plenty of water.
- · Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- · Reference to other sections
- See Section 7 for information on safe handling.
- See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.
- · Protective Action Criteria for Chemicals

· PAC-1:

Substance is not listed.

(Contd. on page 3)

US

Reviewed on 12/03/2015

(Contd. of page 1)



*Printing date 02/16/2017* 

Reviewed on 12/03/2015

(Contd. of page 2)

#### Trade name: DNase/RNase Free Water

#### · PAC-2:

Substance is not listed.

· PAC-3:

Substance is not listed.

#### 7 Handling and storage

- · Handling:
- · Precautions for safe handling No special measures required.
- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- · Storage:
- · Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: None.
- · Specific end use(s) Laboratory reagent

## 8 Exposure controls/personal protection

• Additional information about design of technical systems: No further data; see item 7.

- · Control parameters
- · Components with limit values that require monitoring at the workplace: Not required.
- · Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment:
- · General protective and hygienic measures:
- The usual precautionary measures for handling chemicals should be followed.
- · Breathing equipment: Not required.
- Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation  $\cdot$  *Material of gloves* 

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

• *Eye protection:* Goggles recommended during refilling.

(Contd. on page 4)

Page 3/8



Printing date 02/16/2017

Trade name: DNase/RNase Free Water

Reviewed on 12/03/2015

(Contd. of page 3)

Physical and chemical proper	ties
· Information on basic physical and o	chemical properties
· General Information	I I I
· Appearance:	
Form:	Liquid
Color:	Clear
· Odor:	Odorless
· Odor threshold:	Not determined.
· pH-value:	Not determined.
· Change in condition	
Melting point/Melting range:	0 °C (32 °F)
Boiling point/Boiling range:	100 °C (212 °F)
· Flash point:	Not applicable.
· Flammability (solid, gaseous):	Not applicable.
· Ignition temperature:	
Decomposition temperature:	Not determined.
• Auto igniting:	Not determined.
Danger of explosion:	Product does not present an explosion hazard.
· Explosion limits:	
Lower:	Not Applicable
Upper:	Not Applicable
· Vapor pressure at 20 °C (68 °F):	23 hPa (17 mm Hg)
• Density at 20 •C (68 •F):	1 g/cm <sup>3</sup> (8.345 lbs/gal)
· Relative density	Not determined.
· Vapor density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with	
Water:	Fully miscible.
· Partition coefficient (n-octanol/wate	er): Not determined.
· Viscosity:	
<i>Dynamic at 20 •C (68 •F):</i>	0.952 mPas
Kinematic:	Not determined.
· Solvent content:	
Organic solvents:	0.0 %
Water:	100.0 %
VOC content:	0.0 g/l / 0.00 lb/gl
$\cdot$ Other information	No further relevant information available.

(Contd. on page 5)

US

## Page 4/8



*Printing date 02/16/2017* 

Reviewed on 12/03/2015

Trade name: DNase/RNase Free Water

(Contd. of page 4)

Page 5/8

## 10 Stability and reactivity

- · *Reactivity* No further relevant information available.
- · Chemical stability
- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

## 11 Toxicological information

- · Information on toxicological effects
- Acute toxicity:
- · Primary irritant effect:
- on the skin: No irritant effect.
- on the eye: No irritating effect.
- Sensitization: No sensitizing effects known.
- · Additional toxicological information:

When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

The substance is not subject to classification.

#### · Carcinogenic categories

· IARC (International Agency for Research on Cancer)

Substance is not listed.

· NTP (National Toxicology Program)

Substance is not listed.

· OSHA-Ca (Occupational Safety & Health Administration)

Substance is not listed.

## 12 Ecological information

- · Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes: Generally not hazardous for water
- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- *vPvB*: Not applicable.

(Contd. on page 6)

US



*Printing date 02/16/2017* 

Trade name: DNase/RNase Free Water

· Other adverse effects No further relevant information available.

## 13 Disposal considerations

- · Waste treatment methods
- · *Recommendation:* Smaller quantities can be disposed of with household waste.
- · Uncleaned packagings:
- · Recommendation:
- Dispose of container in acoordance with local/regional/national and international recommendations.
- · Recommended cleansing agent: Water, if necessary with cleansing agents.

UN-Number DOT, ADN, IMDG, IATA	not regulated	
UN proper shipping name DOT, ADN, IMDG, IATA	not regulated	
Transport hazard class(es)		
DOT, ADN, IMDG, IATA		
Class	not regulated	
Packing group		
DOT, IMDG, IATA	not regulated	
Environmental hazards:	Not applicable.	
Special precautions for user	Not applicable.	
Transport in bulk according to Annex II	I of	
MARPOL73/78 and the IBC Code	Not applicable.	

# 15 Regulatory information

· Safety, health and environmental regulations/legislation specific for the substance or mixture · Sara

· Section 355 (extremely hazardous substances):

Substance is not listed.

· Section 313 (Specific toxic chemical listings):

Substance is not listed.

• TSCA (Toxic Substances Control Act):

Substance is listed.

Page 6/8

Reviewed on 12/03/2015

(Contd. of page 5)

(Contd. on page 7)

US



Printing date 02/16/2017

Reviewed on 12/03/2015

Trade name: DNase/RNase Free Water

(Contd. of page 6)

Page 7/8

Proposition 65
 Chemicals known to cause cancer:

Substance is not listed.

· Chemicals known to cause reproductive toxicity for females:

Substance is not listed.

· Chemicals known to cause reproductive toxicity for males:

Substance is not listed.

· Chemicals known to cause developmental toxicity:

Substance is not listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

Substance is not listed.

· TLV (Threshold Limit Value established by ACGIH)

Substance is not listed.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

Substance is not listed.

· GHS label elements Void

· Hazard pictograms Void

· Signal word Void

· Hazard statements Void

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS:

Zymo Research Corp. Safety Department 17062 Murphy Ave. Irvine, CA 92614 USA Phone: 1-949-679-1190 or 1-888-882-9682 · Contact: sds@zymoresearch.com · Date of preparation / last revision 02/16/2017 / -· Abbreviations and acronyms: ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road) IMDG: International Maritime Code for Dangerous Goods DOT: US Department of Transportation IATA: International Air Transport Association ACGIH: American Conference of Governmental Industrial Hygienists EINECS: European Inventory of Existing Commercial Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society) NFPA: National Fire Protection Association (USA)

(Contd. on page 8)

US



Printing date 02/16/2017

Reviewed on 12/03/2015

Trade name: DNase/RNase Free Water

(Contd. of page 7)

Page 8/8

HMIS: Hazardous Materials Identification System (USA) VOC: Volatile Organic Compounds (USA, EU) PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative NIOSH: National Institute for Occupational Safety OSHA: Occupational Safety & Health TLV: Threshold Limit Value PEL: Permissible Exposure Limit REL: Recommended Exposure Limit