

5171 Wilfong Road Memphis, TN 38134 USA

Telephone: 901-382-8716 Fax: 901-333-8223

Email: info@meridianlifescience.com www.MeridianLifeScience.com

## CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial content	Important Note:	Centrifuge before	opening to ensur	e complete recov	ery of vial contents
---	-----------------	-------------------	------------------	------------------	----------------------

**Catalog #:** C01936M **Lot #:** n/a

**Description:** MAb to E. histolytica & E. dispar

Monoclonal Antibody to Entamoeba histolytica and Entamoeba dispar

**Specificity:** Recognizes *Entamoeba histolytica* and *Entamoeba dispar* trophozoites.

Host Animal: Mouse Isotype: IgG2a

**Source:** Ascites

**Immunogen:** Purified *E. histolytica* trophozoites.

Format: Purified, Liquid

**Purification:** > 90% pure. Protein A Chromatography

**Concentration:** (Lot Specific) mg/mL (OD280nm,  $E^{0.1\%} = 1.3$ )

**Buffer:** 0.01 M Phosphate Buffered Saline, pH 7.2

This product contains no stabilizing proteins.

**Preservative:** 0.1% Sodium Azide

Signature

**Applications:** Suitable for use in ELISA and Lateral Flow Assay. Each laboratory should determine an optimum working

titer for use in its particular application. Other applications have not been tested but use in such assays

should not necessarily be excluded.

Capture Detection

C01936M C01936M (Self-Pairs)

**Storage:** Upon receipt, store at 2-8°C. **DO NOT FREEZE**.

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

21 July 2017

Date