

5171 Wilfong Road Memphis, TN 38134 USA

Telephone: 901-382-8716 Fax: 901-333-8223

Email: info@meridianlifescience.com www.MeridianLifeScience.com

## CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** H86177M **Lot #:** 3C06716

**Description:** MAb to Serum Amyloid A

Monoclonal Antibody to Human Serum Amyloid A (SAA)

**Specificity:** Recognizes native and recombinant Human SAA.

**Host Animal:** Mouse Isotype:  $IgG_1$ 

**Source:** Ascites

Immunogen: Human recombinant SAA, conjugated with carrier protein.

Format: Purified, Liquid

**Purification:** > 90% pure (SDS-PAGE). Protein A Chromatography

**Concentration:** 5.5 mg/mL (OD280nm,  $E^{0.1\%} = 1.4$ )

**Buffer:** Phosphate Buffered Saline, pH 7.4

**Preservative:** 0.09% Sodium Azide

**Applications:** Suitable for use in ELISA and Western Blot. Each laboratory should determine an optimum working titer for

use in its particular application. Other applications have not been tested but use in such assays should not

necessarily be excluded.

Recommended pairs for sandwich immunoassay:

 Capture
 Detection

 H86177M
 H86180M

 H86180M
 H86177M

 H86177M
 H86178M

 H86178M
 H86177M

**Storage:** Store at 2-8°C.

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

13 September 2016

Signature Date