



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: K8A069M **Lot #:** 6C08812

Description: MAb to CD37
Monoclonal Antibody to CD37

Specificity: CD37 antigen (human), 40-52 kDa glycoprotein.

Clone: IPO-24

Host Animal: Mouse **Isotype:** IgG_{2b}

Source: Ascites

Immunogen: Produced against a B cell line.

Format: Purified, Liquid

Purification: Protein A Chromatography

Concentration: 1 mg/mL

Buffer: 1X PBS, pH 7.2

Preservative: None

Application: Suitable for use in Immunoprecipitation, ELISA (1:1,000), IHC (paraffin sections and frozen sections). Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Short term (up to 2 months) store at 2–8°C. Long term, aliquot and store at -80°C. Avoid multiple freeze/thaw cycles.

Signature

17 Dec 2015

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY