



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com
www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Catalog #: L66644G **Lot #:** 8K30717

Description: Goat anti Human C4
Goat Antibody to Human Complement Component 4
Turbidimetric Grade

Specificity: Monospecific by Immunoelectrophoresis (IEP) against normal human plasma and 2X concentrated normal human serum.

Host Animal: Goat

Immunogen: Purified Human C4

Format: Monospecific Antiserum, Liquid

Purification: Delipidation followed by solid phase absorption.
Product is 0.2 µm filtered.

Concentration: Total Protein: 6.6 g/dL
Becker Titer: 7.4 mg/mL

Buffer: 0.05 M Tris, 0.5 M Sodium Chloride, pH 7.5.

Preservative: 0.09% Sodium Azide

Applications: Suitable for use in Immunoelectrophoresis and Turbidimetry. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store at 2-8°C. Long term store at -20°C.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Diana Dun

03Nov2017

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY