



## Normal Donkey Serum

Art. no PSB-T0Q10H-10  
Size 10 mL

<b>Source</b>	Adult donkeys
<b>Purification</b>	Processed to minimize hemolysis and 0.2 µm filtered
<b>Applications</b>	Block
<b>Format</b>	Serum
<b>Storage &amp; Handling</b>	The serum is supplied as 10.0 mL with 0.1% NaN <sub>3</sub> as a preservative. Store at 2-8°C or aliquot and store at or below -20°C for long term storage. Reagent is stable for the period shown on the label if stored as directed.
<b>Warning</b>	Reagent contains sodium azide. Please refer to product specific M)SDS.

For troubleshooting or questions regarding the product, please contact our support team at [support@nordicbiosite.com](mailto:support@nordicbiosite.com)

For RESEARCH USE ONLY [RUO]