



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: Q2B101G **Lot #:** 24L34101

Description: Goat A γ Amyloid Pre.Pro. Frame (C-Term)
Goat Antibody to Human Amyloid Precursor Protein Frameshift Mutant C-Terminal

Specificity: Amino acid sequence 339-348 of human amyloid precursor protein frameshift mutant.

Host Animal: Goat

Immunogen: Synthetic Peptide: ${}_{339}\text{RGRTSSKELA}_{348}$

Format: Neat, Liquid

Purification: Not Applicable

Concentration: Not Determined

Buffer: Not Determined

Preservative: 0.1% Sodium Azide

Applications: ELISA using immunizing peptide ($\geq 1:5,000$).
IHC using formalin-fixed, paraffin-embedded sections ($\geq 1:300$).
Western Blot ($\geq 1:3,000$).
Suitability for use in IHC using paraformaldehyde-fixed, frozen sections has not been determined but use in such assays should not necessarily be excluded.

Storage: Short term store at 2-8°C. Long term store at -20°C/

Warnings: This product contains sodium azide, which has been classified as Xn (Harmful) in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.

Signature

20 September 2013

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY