



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** R01523 **Lot #:** 1F16517

**Description:** *B. burgdorferi* VlsE, Recombinant  
*Borrelia burgdorferi* Variable major protein-like sequence expressed (VlsE), Recombinant.

Recombinant *Borrelia burgdorferi* VlsE antigen produced from the complete ORF of the variable major protein-like sequence expressed. Contains a His-tag at the N-terminus. Relative molecular weight calculated from the amino acid sequence is 39,926.9 Da.

**Source:** *E. coli*

**Format:** Purified, Liquid

**Purification:** 95% pure (SDS-PAGE).

**Concentration:** 1.13 mg/mL (Bradford Method)

**Buffer:** 20 mM Phosphate, pH 7, 1 M Sodium Chloride, 0.1% Polyoxyethylene(10)Tridecylether.

**Preservative:** None

**Applications:** Suitable for IgG and IgM detection in ELISA (0.32 µg/mL), Dot Blot and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage:** Store at -20°C. Aliquot to avoid multiple freeze/thaw cycles.

**Inactivation:** Not Applicable

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Corry Thomas

14 Jun 2017

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**