



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: R01764 **Lot #:** n/a

Description: Enterovirus A EV71 Recomb.
Enterovirus A EV71 Recombinant

Product contains a GST fusion partner.
Molecular Weight: 60 kDa

Source: *E. coli*

Format: Purified, Liquid

Purification: > 85% Purity. Multi-step procedure including affinity chromatography.

Concentration: (Lot Specific) mg/mL (OD280nm, E^{0.1%} = 1.40)

Buffer: Tris Buffer, pH 8.0

Preservative: 0.09% Sodium Azide

Applications: Suitable for use as a control in Lateral Flow Assays and as an immunogen. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store at 2-8°C for < 1 week. For long term storage, aliquot and store at -20°C to avoid multiple freeze/thaw cycles.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Signature

26 January 2018

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY