



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** R86350 **Lot #:** 1A01718

**Description:** Influenza B (B/Victoria/504/00)

**Source:** Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza B virus.  
Strain: B/Victoria/504/00

**Format:** Purified, Liquid

**Purification:** > 90% pure (SDS-PAGE). Ultracentrifugation using 10-40% sucrose gradient.

**Concentration:** 0.9 mg/mL

**Buffer:** 0.05 M Tris-HCl, pH 8.0 containing 0.1 M Sodium Chloride, 5 mM EDTA.

**Preservative:** 0.09% Sodium Azide and 0.005% Thimerosal.

**Applications:** Serological studies of influenza B virus, immunogen for antibody production. Has been tested with anti-influenza B monoclonal antibodies in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage:** Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.

**Inactivation:** 0.005% Thimerosal and beta propiolactone.

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

**References:** The reference listed below is for research purposes only:  
Yan, X., et al., (2005), "Multiplexed Flow Cytometric Immunoassay for Influenza Virus Detection and Differentiation", Anal. Chem., **77**: 7673-7678.

*Brenda Dinn*

18Jan2018

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**