



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road

Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: Y01235B **Lot #:** 8K32415

Description: Buprenorphine-BSA Conjugate
Buprenorphine-Bovine Serum Albumin (BUP-BSA) Conjugate
Bovine Serum Albumin Conjugated

Format: BSA, Liquid

Concentration: 8.0 mg/mL

Buffer: PBS, pH 7.2

Preservative: 0.02% Sodium Azide

Applications: Specific methodologies have not been tested using this product.

Storage: Short term, store at 2–8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Signature

19 May 2016

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

001