



Meridian
Life Science, Inc.
Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA
Telephone: 901-382-8716
Fax: 901-333-8223
Email: info@meridianlifescience.com
www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: Y01285B **Lot #:** 26L34215

Description: Barbiturate BSA Conjugate
Barbiturate Bovine Serum Albumin Conjugate

Source: Not Applicable

Format: BSA, Liquid

Concentration: 16.7 mg/mL

Buffer: 10 mM PBS, pH 7.4

Preservative: 0.1% Sodium Azide

Applications: Suitable for use in ELISA and Lateral Flow assays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended antibody-antigen pair for competitive inhibition immunoassays:

Antigen-Capture

Antibody-Labeling

Y01285B

G01285M

Y01285B

G01283M

Storage: Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Robert Ott

Signature

23 May 2016

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY