

5171 Wilfong Road Memphis, TN 38134

Telephone: 901-382-8716 Fax: 901-333-8223

Email: info@meridianlifescience.com www.MeridianLifeScience.com

## CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** Y01297B **Lot #:** 7L36117

**Description:** Morphine BSA Conjugate

Morphine Bovine Serum Albumin (BSA) Conjugate

**Source:** Not Applicable

Format: BSA, Liquid (Hapten/Protein Ratio ~70:1; Link Position 3)

**Concentration:** 8.95 mg/mL

**Buffer:** 10 mM Phosphate Buffered Saline, pH 7.4

**Preservative:** 0.1% Sodium Azide

**Applications:** Suitable for use in ELISA and Lateral Flow assays. Each laboratory should determine an optimum working

titer for use in its particular application. Other applications have not been tested but use in such assays

should not necessarily be excluded.

Recommended antigen-antibody pair for competitive inhibition immunoassays:

Antigen-Capture Antibody-Labeling

Y01297B G01298M Y01297B G01299M Y01297B G01300M

**Storage:** Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

**Inactivation:** Not Applicable

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Brenda Dum

28Dec2017

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY