



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
Telephone: 901-382-8716
Fax: 901-333-8223
Email: info@meridianlifescience.com
www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: Y82135P **Lot #:** 6C06005

Description: Methamphetamine-S-HRP Conjugate
Methamphetamine-S Horseradish Peroxidase (HRP) Conjugate
Horseradish Peroxidase conjugated

Format: HRP, Liquid

Purification: Not applicable.

Concentration: Not determined.

Buffer: Tris, pH 7.6 with protein stabilizers

Preservative: 0.005% Thimerosal and other additives
Note: Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.

Applications: Specific methodologies have not been tested using this product.

Storage: Store at 2-8°C.

Signature

27 May 2016

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY