

5171 Wilfong Road Memphis, TN 38134 USA

Telephone: 901-382-8716 Fax: 901-333-8223

Email: info@meridianlifescience.com www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #:	Y82290P	Lot #:	10B04606

Description: Cotinine-3-HRP Conjugate
Horseradish Peroxidase (HR

Horseradish Peroxidase (HRP) Conjugate

Source: Not Applicable

Important Note:

Format: HRP, Liquid

Purification: HRP coupled to Cotinine at position 3.

Concentration: Not Determined

Buffer: 0.05 M Tris containing protein stabilizers.

Preservative: 0.005% Thimerosal and other additives.

Warning: Use of sodium azide as a preservative will substantially inhibit the enzyme activity of

horseradish peroxidase.

Applications: Specific methodologies have not been tested using this product.

Storage: Store at 2-8°C.

Robert OH

27 May 2016