



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** Y82451P **Lot #:** 4B04916

**Description:** Cotinine-4-HRP Conjugate  
Horseradish Peroxidase (HRP) Conjugate

**Source:** Not Applicable

**Format:** HRP, Liquid

**Purification:** HRP coupled to Cotinine at position 4.

**Concentration:** Not Determined

**Buffer:** 0.05 M Tris containing protein stabilizers.

**Preservative:** 0.005% Thimerosal and other additives.  
**Warning:** Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.

**Applications:** Specific methodologies have not been tested using this product.

**Storage:** Store at 2–8°C. **DO NOT FREEZE!**

Signature

27 May 2016

Date

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**